

***PROGRAMMING GUIDE***

**BÓSA**





## Table of contents

Table of contents .....	3
0 Introduction .....	5
1 Using the Ethernet Interface .....	7
2 Using the GPIB Interface .....	9
3 Using the Macro Editor Tool .....	11
3.1 Macro Editor Command Summary .....	11
3.2 Macro Editor Command Definitions .....	15
3.2.1 Array Functions .....	16
3.2.2 Conditional Statements and functions .....	17
3.2.3 BOSA control instruments .....	20
3.2.4 Mathematical functions .....	21
3.2.5 Date and Time functions .....	27
3.2.6 String functions .....	30
3.2.7 Explicit conversion functions .....	31
4 BOSA SCPI Commands .....	33
4.1 BOSA SCPI Command Summary .....	33
4.2 BOSA SCPI Command Definitions .....	42
4.2.1 INSTRument Subsystem Commands .....	42
4.2.2 DISPlay Subsystem Commands .....	43
4.2.3 SENSE Subsystem Commands .....	46
4.2.4 INPut Subsystem Commands .....	54
4.2.5 CALCulate Subsystem Commands .....	56
4.2.6 TRACe Subsystem Commands .....	64
4.2.7 FORMat Subsystem Commands .....	65
4.2.8 MMEMory Subsystem Commands .....	66
5 Programming examples .....	69
5.1 Ethernet Interface .....	69
5.2 Macro Editor Tool .....	70
5.2.1 Performing auto measurements of several peaks using markers ..	70
5.2.2 Wavelength shift measurement on a CW laser .....	72
5.2.3 Controlling an external device via GPIB .....	73



# 0 Introduction

This programming guide gives information on how to control and program the BOSA equipment. The user is allowed to programmatically control the BOSA in three different ways:

- The Macro Editor Tool: it is a built-in functionality in the BOSA graphic user interface. It allows the user to program automatic measurements, set analysis functions and save results. Through simple scripts all BOSA functionalities are accessible. The Macro Editor Tool also enables communication with any external equipment via GPIB acting the BOSA as controller. The specific language syntax and SCPI commands are described in sections 3 and 4 in this manual.
- GPIB: the General Purpose Interface Bus is a common used way of controlling laboratory equipment. The GPIB interface of the BOSA allows remote control of the system allowing to automate measurements and to integrate the BOSA in your measurement set up. GPIB configuration and SCPI commands are described in sections 2 and 4 in this manual.
- Ethernet: the Ethernet interface is a standard interface used for personal computers to communicate over the network. The Ethernet connector on the rear panel of the BOSA allows remote communication through standard network cable. Network parameters and SCPI commands are described in sections 1 and 4 in this manual. Default password is 1234.

Macro Editor Tool, GPIB and Ethernet Interface use the same set of SCPI commands to operate the BOSA. A detailed description of the commands can be found in section 3.



# 1 Using the Ethernet Interface

When using the Ethernet interface the BOSA acts as a server. To operate, connect the Ethernet connector on the back of the BOSA to the LAN-network by using a standard network cable. The IP address of the BOSA should be configured to DHCP or manual (please refer to section 5.2 in the BOSA User Guide for details).

The port is 10000 on the TCP/IP protocol. Any command sent has to be followed by a read command. For non query commands the answer is "OK" or an error message. The error messages can be one of the followings: command error, parameter error or unit error.

The BOSA server will only admit one client at a time. Other clients trying to connect to the BOSA will not be taken into account.





## 2 Using the GPIB Interface

The GPIB (General Purpose Interface Bus) is an interface used for communication between personal computers and laboratory instruments, such as the BOSA. The IEEE-488 connector on the back of the BOSA allows for remote operation through a standard GPIB cable. The connector is a standard female 24-pin IEEE-488 connector to use with a standard shielded IEEE-488 cable.

Before you can operate the BOSA through GPIB interface you must assign the BOSA a device address that is unique from the other equipment attached to the bus.

### Setting the Device Address

Under the Settings menu, select GPIB. The GPIB Configuration window appears. Select there a GPIB primary address for the BOSA different from the controller computer or the other instruments' addresses.

The default GPIB address of the BOSA is 4.

### Returning to Local Control

When the BOSA is in remote control, the text of the Stop/Run button in the main screen of the BOSA will change to Local and all the main interface of the BOSA will be disabled. Press Local button if you wish to return the instrument to local control.

The complete set of commands that can be sent through GPIB interface are listed in section 3 of this manual. Additionally to that set of commands the BOSA will respond to the standard Identification query.

Identification query:

<b>Syntax:</b>	<b>*IDN?</b>
<b>Description:</b>	Returns the system identification string of the BOSA.
<b>Example of response:</b>	ARAGON-PHOTONICS,BOSA-C,AC122201151010,V1.3.42 (Manufacturer: Aragon photonics, Model: Bosa-c, Hardware version: AC122201151010, Firmware version: V1.3.42)

Operation complete query:

**Syntax:** `*OPC?`

**Description:** Returns an ASCII "1" when all pending overlapped operations have been completed.

**Example of response:** 1

## 3 Using the Macro Editor Tool

The *Macro Editor Tool* is a programming function that allows the user to program BOSA measurements conditions, set analysis functions and save results. It also enables communications with any external equipment acting the BOSA as controller. In the Macro editor are defined the following specific command categories:

- Variables and Arrays
- Conditional Statements and Functions
- BOSA Control Functions
- Mathematical Functions
- Date and Time Functions
- String Functions
- Explicit Conversion Functions

### 3.1 Macro Editor Command Summary

The following table gives an overview of the commands used in the Macro Editor, a description and a page reference for more detailed information about the particular command.

ARRAY FUNCTIONS		
Array <b>array</b> (int <i>length</i> )	Creates a new array with the specified <i>length</i> .	16
Void <b>sort</b> (string <i>arrayName</i> )	Sort the given array, identified by <i>arrayName</i> , into ascending order.	16
Void <b>reverse</b> (string <i>arrayName</i> )	Reverses the items of the given array, identified by <i>arrayName</i> .	16
Int <b>length</b> (string <i>arrayName</i> )	Returns the length of the given array, identified by <i>arrayName</i> .	16
Void <b>trim</b> (string <i>arrayName</i> )	Reduces the size of the given array, identified by <i>arrayName</i> , by eliminating the null values.	17

Void <b>plot</b> (string arrayName)	Plots the value of the arrayName in a chart	17
Void <b>plotXY</b> (string arrayNameX,string arrayNameY)	Plots the pair of values given by arrayNameX and arrayNameY in a XY chart.	17
<b>CONDITIONAL STATEMENTS AND FUNCTIONS</b>		
<b>if</b> (bool <i>condition</i> ) ... <b>endif</b>	Performs a conditional statement.	17
double <b>iif</b> (bool <i>condition</i> , double <i>val1</i> , double <i>val2</i> )	Returns <i>val1</i> if <i>condition</i> is true, otherwise returns <i>val2</i> .	18
double <b>case</b> (bool <i>condition1</i> , double <i>val1</i> , bool <i>condition2</i> , double <i>val2</i> )	Returns the value after the first condition evaluated as true.	18
<b>repeat</b> ... <b>while</b> (bool <i>condition</i> )	If <i>condition</i> is true, the code execution jumps to the line after the repeat statement.	18
<b>break</b>	Jumps out of the actual repeat-while loop.	19
Void <b>wait</b> (int <i>value</i> )	Forces the BOSA to wait for the specified time before executing the next command. <i>Value</i> is given in milliseconds.	19
<b>BOSA CONTROL INSTRUMENTS</b>		
String <b>scpi</b> (string <i>command</i> )	Sends a <i>scpi command</i> . If it is a query command the answer is returned, otherwise it returns OK. An error string is returned if the command, parameters or units are not correct.	20
Void <b>sendgpib</b> (int <i>address</i> , string <i>command</i> )	Sends a GPIB <i>command</i> to an external instrument configured at the specified primary <i>address</i> . The command is specific for the external instrument.	20
String <b>querygpib</b> (int <i>address</i> , string <i>command</i> )	Sends a query GPIB <i>command</i> to an external instrument configured at the specified primary <i>address</i> . The command is specific for the external instrument.	20
Void <b>fwrite</b> (string <i>file</i> , string <i>value</i> )	Writes the <i>value</i> in the desired <i>file</i> .	21
Void <b>fwriteln</b> (string <i>file</i> , string <i>value</i> )	Writes a line with the <i>value</i> in the desired <i>file</i> .	21
Void <b>setvar</b> (string <i>name</i> , object <i>value</i> )	Creates a new variable identified by <i>name</i> , and initiates its <i>value</i> .	21
<b>MATHEMATICAL FUNCTIONS</b>		
Double <b>e</b> ()	Returns the value of the e constant.	21

Double <b>pi</b> ()	Returns the value of the Pi constant.	21
Double <b>c</b> ()	Returns the speed of light in vacuum.	22
Double <b>sin</b> (double <i>value</i> )	Returns the sine of the specified <i>value</i> .	22
Double <b>cos</b> (double <i>value</i> )	Returns the cosine of the specified <i>value</i> .	22
Double <b>tan</b> (double <i>value</i> )	Returns the tangent of the specified <i>value</i> .	22
Double <b>sec</b> (double <i>value</i> )	Returns the secant of the specified <i>value</i> .	23
Double <b>csc</b> (double <i>value</i> )	Returns the cosecant of the specified <i>value</i> .	23
Double <b>cot</b> (double <i>value</i> )	Returns the cotangent of the specified <i>value</i> .	23
Double <b>asin</b> (double <i>value</i> )	Returns the arcsine of the specified <i>value</i> .	23
Double <b>acos</b> (double <i>value</i> )	Returns the arccosine of the specified <i>value</i> .	23
Double <b>atan</b> (double <i>value</i> )	Returns the arctangent of the specified <i>value</i> .	24
Double <b>sinh</b> (double <i>value</i> )	Returns the hyperbolic sine of the specified <i>value</i> .	24
Double <b>cosh</b> (double <i>value</i> )	Returns the hyperbolic cosine of the specified <i>value</i> .	24
Double <b>tanh</b> (double <i>value</i> )	Returns the hyperbolic tangent of the specified <i>value</i> .	24
Double <b>abs</b> (double <i>value</i> )	Returns the absolute value of the specified <i>value</i> .	24
Double <b>sqrt</b> (double <i>value</i> )	Returns the square root of the specified <i>value</i> .	25
Int <b>ceil</b> (double <i>value</i> )	Returns the smallest integer greater than the specified <i>value</i> .	25
Int <b>floor</b> (double <i>value</i> )	Returns the largest integer smaller than the specified <i>value</i> .	25
Double <b>exp</b> (double <i>value</i> )	Returns e raised to a power indicated by <i>value</i> .	25
Double <b>log10</b> (double <i>value</i> )	Returns the decimal logarithmic of the specified <i>value</i> .	25
Double <b>log</b> (double <i>value1</i> , [double <i>base</i> ])	Returns the logarithmic of the specified value, in the specified <i>base</i> .	26
Double <b>max</b> (double <i>value1</i> , double <i>value2</i> )	Returns the largest value of the two specified values.	26

Double <b>min</b> (double <i>value1</i> , double <i>value2</i> )	Returns the smallest value of the two specified values.	26
Double <b>pow</b> (double <i>value1</i> , double <i>value2</i> )	Returns <i>value1</i> raised to <i>value2</i> .	26
Double <b>round</b> (double <i>value1</i> , int <i>value2</i> )	Returns <i>value1</i> rounded with the number of decimal digits specified by <i>value2</i> .	27
<b>DATE &amp; TIME FUNCTIONS</b>		
DateTime <b>now</b> ()	Returns the actual date time.	27
DateTime <b>today</b> ()	Returns the actual date.	27
DateTime <b>mindate</b> ()	Returns the minimum value for the date (00:00:00.0000000, 1 <sup>st</sup> of January, 0001).	27
DateTime <b>maxdate</b> ()	Returns the maximum value for the date (23:59:59.9999999, 31 <sup>st</sup> of December, 9999).	28
String <b>monthname</b> ()	Returns the name of the actual month.	28
DateTime <b>adddays</b> (dateTime <i>date</i> , int <i>days</i> )	Adds the specified amount of <i>days</i> to the specified <i>date</i> .	28
DateTime <b>addmonths</b> (dateTime <i>date</i> , int <i>months</i> )	Adds the specified amount of <i>months</i> to the specified <i>date</i> .	28
DateTime <b>addyears</b> (dateTime <i>date</i> , int <i>years</i> )	Adds the specified amount of <i>years</i> to the specified <i>date</i> .	28
DateTime <b>addhours</b> (dateTime <i>time</i> , int <i>hours</i> )	Adds the specified amount of <i>hours</i> to the specified <i>time</i> .	28
DateTime <b>addminutes</b> (dateTime <i>time</i> , int <i>minutes</i> )	Adds the specified amount of <i>minutes</i> to the specified <i>time</i> .	29
DateTime <b>addseconds</b> (dateTime <i>time</i> , int <i>seconds</i> )	Adds the specified amount of <i>seconds</i> to the specified <i>time</i> .	29
<b>STRING FUNCTIONS</b>		
String <b>concat</b> (string <i>value1</i> , string <i>value2</i> , string <i>value3</i> ...)	Concatenates the strings <i>value1</i> , <i>value2</i> ....	30
String <b>formatnum</b> (double <i>value</i> , string <i>format</i> )	Returns a string with the <i>value</i> specified in the <i>format</i> specified.	30
String <b>formatdate</b> (dateTime <i>value</i> , string <i>format</i> )	Returns the date Time <i>value</i> specified in the <i>format</i> specified.	30
<b>EXPLICIT CONVERSION FUNCTIONS</b>		
Double <b>todouble</b> (object <i>value</i> )	Converts the specified <i>value</i> to a double.	31
Int <b>toint</b> (object <i>value</i> )	Converts the specified <i>value</i> to an integer.	31

Long <b>tolong</b> (object <i>value</i> )	Converts the specified <i>value</i> to a long.	31
Unsigned int <b>touint</b> (object <i>value</i> )	Converts the specified <i>value</i> to an unsigned integer.	31
Unsigned long <b>toulong</b> (object <i>value</i> )	Converts the specified <i>value</i> to an unsigned long.	31
DateTime <b>todatetime</b> (object <i>value</i> )	Converts the specified <i>value</i> to a date Time.	32
String <b>tostring</b> (object <i>value</i> )	Converts the specified <i>value</i> to a string.	32

## 3.2 Macro Editor Command Definitions

A new variable is declared just by assigning it a value. There is no need to define the type of the new variable and both functions and operators will perform automatic type conversion if needed.

Operator	Description	Example	Result
+	Addition	4 + 1	5
-	Subtraction	4 - 2	3
*	Multiplication	4 * 3	12
/	Division	4 / 3	1.3333
%	Modulo	4 % 3	1
&	Concatenation	1 & 2	12

Others:

Operator	Description	Example	Result
==	Comparison	4 == 4	True
!=	Negative comparison	4 != 3	True
//	Comment		

Some examples:

- a=1
- b=2

- $c = a + b$  (result:  $c = 3$ )
- $c = a \& b$  (result:  $c = 12$ )
- $d = "3"$
- $f = a/d$  (result:  $f = 0.333$ )

### 3.2.1 Array Functions

#### 3.2.1.1 Array

<b>Syntax</b>	<b>Array array(int <i>length</i>)</b>
<b>Description</b>	Creates a new array with the specified length.
<b>Argument</b>	'length' specifies the size of the array.
<b>Example</b>	myArray = array(10)
<b>Comments</b>	The first element is indexed with 0. When created, all elements are initialized to null values. To access an element: myArray[8]

#### 3.2.1.2 Sort

<b>Syntax</b>	<b>Void sort(arrayName)</b>
<b>Description</b>	Sorts a given array into ascending order.
<b>Argument</b>	'arrayName' is the name of the array to be sorted.
<b>Example</b>	sort(myArray)
<b>Comments</b>	Only numeric arrays can be used.

#### 3.2.1.3 Reverse

<b>Syntax</b>	<b>Void reverse(arrayName)</b>
<b>Description</b>	Reverses the items of a given array.
<b>Argument</b>	'arrayName' is the name of the array to be reversed.
<b>Example</b>	reverse(myArray)
<b>Comments</b>	None.

#### 3.2.1.4 Length

<b>Syntax</b>	<b>Int length(arrayName)</b>
<b>Description</b>	Returns the length of the given array.
<b>Argument</b>	'arrayName' is the name of the array.



<b>Example</b>	Length(myArray)
<b>Comments</b>	None.

### 3.2.1.5 Trim

<b>Syntax</b>	<b>Void trim(arrayName)</b>
<b>Description</b>	Reduces the size of the given array by eliminating the null values.
<b>Argument</b>	'arrayName' is the name of the array.
<b>Example</b>	trim(myArray)
<b>Comments</b>	Note that null value is not the same as zero value or empty (in strings).

### 3.2.1.6 Plot

<b>Syntax</b>	<b>Void plot(arrayName)</b>
<b>Description</b>	Plots the values of the array in a chart
<b>Argument</b>	'arrayName' is the name of the data array to display.
<b>Example</b>	plot(myArray)

### 3.2.1.7 PlotXY

<b>Syntax</b>	<b>Void plotXY(arrayNameX,arrayNameY)</b>
<b>Description</b>	Plots the pair of values of the two arrays in a XY chart
<b>Argument</b>	'arrayNameX' is the name of the data array to be displayed in the X axis and 'arrayNameY' is the data array to be displayed in the Y axis .
<b>Example</b>	plotXY(myArrayX,myArrayY)
<b>Comments</b>	ArrayX and arrayY lengths must be the same

## 3.2.2 Conditional Statements and functions

### 3.2.2.1 If

<b>Syntax</b>	<b>if(bool <i>condition</i>)</b>  <b>...</b> <b>Endif</b>
<b>Description</b>	Performs a conditional statement.

<b>Argument</b>	'condition' can be a boolean value or expression.
<b>Example</b>	<pre> if(bool <i>condition</i>)     <i>statement1</i>     if(<i>condition2</i>)         <i>statement3</i>         ...     endif     <i>statement2</i>     ... Endif </pre>
<b>Comments</b>	Several condition statements can be nested.

### 3.2.2.2 Double iif

<b>Syntax</b>	<b>double iif(bool <i>condition</i>, double <i>val1</i>, double <i>val2</i>)</b>
<b>Description</b>	It returns the object val1 if condition is true, otherwise object val2.
<b>Argument</b>	<p>'condition' can be a boolean value or expression.</p> <p>'val1' is the result value if condition is true.</p> <p>'val2' is the result value if condition is false.</p>
<b>Example</b>	myMajor = iif(a>b,a,b)
<b>Comments</b>	It is expected that, at least one of the conditions is true.

### 3.2.2.3 Double case

<b>Syntax</b>	<b>double case(bool <i>condition1</i>, double <i>val1</i>, bool <i>condition2</i>, double <i>val2</i>)</b>
<b>Description</b>	Returns the parameter after the first condition parameter that evaluates to true.
<b>Argument</b>	<p>'condition1' can be a boolean value or expression.</p> <p>'val1' is the returned value if condition1 is true.</p> <p>'condition2' can be a boolean value or expression.</p> <p>'val2' is the returned value if condition2 is true.</p>
<b>Example</b>	myCase = case(a<0,-1,a=0,0,a>0,1)
<b>Comments</b>	None.

### 3.2.2.4 Repeat-while

<b>Syntax</b>	<pre> repeat ... while(bool <i>condition</i>) </pre>
---------------	--

<b>Description</b>	If condition is true, the code execution jumps to the line after the repeat statement.
<b>Argument</b>	'condition' can be a boolean value or expression.
<b>Example</b>	<pre> repeat     <i>statement1</i>     <i>statement2</i>     if(<i>condition2</i>)         repeat             <i>statement4</i>             ...         while(<i>condition3</i>)     break endif     <i>statement3</i> ... while(bool <i>condition</i>) </pre>
<b>Comments</b>	Each statement is executed at least once.

### 3.2.2.5 Break

<b>Syntax</b>	<b>Break</b>
<b>Description</b>	Jumps out of the actual repeat-while loop
<b>Argument</b>	None.
<b>Example</b>	<pre> repeat     <i>statement 1</i>     if(<i>condition1</i>)         break     endif while(<i>condition2</i>) </pre>
<b>Comments</b>	None.

### 3.2.2.6 Wait

<b>Syntax</b>	<b>Void wait(int <i>value</i>)</b>
<b>Description</b>	Forces the BOSA to wait for the specified time in milliseconds before executing the next command.
<b>Argument</b>	'value' specifies the number of milliseconds to wait.
<b>Example</b>	wait(1000)
<b>Comments</b>	None.

### 3.2.3 BOSA control instruments

#### 3.2.3.1 *String scpi(string command)*

<b>Syntax</b>	<b>String scpi(string <i>command</i>)</b>
<b>Description</b>	Executes a scpi command. If it is a query command the answer is returned, otherwise it returns OK. An error string is returned if the command, parameters or units are not correct.
<b>Argument</b>	'command' is the SCPI command to manage the BOSA.
<b>Example</b>	myValue = scpi("sense:wavelength:center?")
<b>Comments</b>	Refer to section 4, <i>BOSA SCPI Commands</i> .

#### 3.2.3.2 *Void sendgpib(int address, string command)*

<b>Syntax</b>	<b>Void sendgpib(int <i>address</i>, string <i>command</i>)</b>
<b>Description</b>	Sends a GPIB <i>command</i> to an external instrument configured at the specified primary <i>address</i> . The command is specific for the external instrument.
<b>Argument</b>	'address' specifies the GPIB primary address of the device to be controlled from the BOSA. 'command' is the command to be sent to the remote instrument.
<b>Example</b>	sendgpib(8,"state on")
<b>Comments</b>	None.

#### 3.2.3.3 *String querygpib(int address, string command)*

<b>Syntax</b>	<b>String querygpib (int <i>address</i>, string <i>command</i>)</b>
<b>Description</b>	Sends a query GPIB <i>command</i> to an external instrument configured at the specified primary <i>address</i> . The command is specific for the external instrument.
<b>Argument</b>	'address' specifies the GPIB primary address of the device to be controlled from the BOSA. 'command' is the query command to be sent to the remote instrument.
<b>Example</b>	querygpib(8,"*IDN?")
<b>Comments</b>	None.

### 3.2.3.4 Void fwrite(string file, string value)

<b>Syntax</b>	<b>Void fwrite(string <i>file</i>, string <i>value</i>)</b>
<b>Description</b>	Writes the value in the desired file.
<b>Argument</b>	'file' specifies the complete path of the file. 'value' contents the string to be wrote in the file.
<b>Example</b>	Fwrite("D:\User Files\BOSA Data Files\myFile.txt","hello")
<b>Comments</b>	If the file does not exist, it will be created. If omitted, default path is 'D:\User Files\'. If the file already exists, lines will be added to the file.

### 3.2.3.5 Void fwriteln(string file, string value)

<b>Syntax</b>	<b>Void fwriteln(string <i>file</i>, string <i>value</i>)</b>
<b>Description</b>	Writes a line with the value in the desired file.
<b>Argument</b>	'file' specifies the complete path of the file. 'value' contents the line to be wrote in the file.
<b>Example</b>	fwriteln("D:\User Files\BOSA Data Files\myFile.txt","hello")
<b>Comments</b>	If the file does not exist, it will be created. If omitted, default path is 'D:\User Files\'. If the file already exists, lines will be added to the file.

### 3.2.3.6 Void setvar(string name, object value)

<b>Syntax</b>	<b>Void setvar(string <i>name</i>, object <i>value</i>)</b>
<b>Description</b>	Creates a new variable, and initiates its value.
<b>Argument</b>	'name 'specifies the name of the new variable. 'value' sets the value for the new variable.
<b>Example</b>	Setvar("myVar",155)
<b>Comments</b>	The new variable can be a double, an integer, a string, etc. Note the difference between: 'setvar("myVar",155)' and 'setvar("myVar","155")'. In the first case, myVar is an integer; in the second case myVar is a string.

## 3.2.4 Mathematical functions

### 3.2.4.1 Double e()

<b>Syntax</b>	<b>Double e()</b>
<b>Description</b>	Returns the value of the e constant.

<b>Argument</b>	None.
<b>Example</b>	<code>myEE = e()*e()</code>
<b>Comments</b>	None.

#### 3.2.4.2 Double pi()

<b>Syntax</b>	<b>Double pi()</b>
<b>Description</b>	Returns the value of the Pi constant.
<b>Argument</b>	None.
<b>Example</b>	<code>myArea = pi()*myRad*myRad</code>
<b>Comments</b>	None.

#### 3.2.4.3 Doble c()

<b>Syntax</b>	<b>Doble c()</b>
<b>Description</b>	Returns the speed of light in vacuum.
<b>Argument</b>	None.
<b>Example</b>	<code>energy = mass*c()*c()</code>
<b>Comments</b>	In SI units: metres/second, m/s.

#### 3.2.4.4 Double sin(double value)

<b>Syntax</b>	<b>Double sin(double value)</b>
<b>Description</b>	Returns the sine of the specified angle.
<b>Argument</b>	'value' sets the angle, in radians.
<b>Example</b>	<code>mySine = sin(myAngle)</code>
<b>Comments</b>	None.

#### 3.2.4.5 Double cos(double value)

<b>Syntax</b>	<b>Double cos(double value)</b>
<b>Description</b>	Returns the cosine of the specified angle.
<b>Argument</b>	'value' sets the angle, in radians.
<b>Example</b>	<code>myCosine = cos(myAngle)</code>
<b>Comments</b>	None.

#### 3.2.4.6 Double tan(double value)

<b>Syntax</b>	<b>Double tan(double value)</b>
<b>Description</b>	Returns the tangent of the specified angle.

<b>Argument</b>	'value' sets the angle, in radians.
<b>Example</b>	myTangent = tan(myAngle)
<b>Comments</b>	None.

#### 3.2.4.7 Double sec(double value)

<b>Syntax</b>	<b>Double sec(double value)</b>
<b>Description</b>	Returns the secant of the specified angle.
<b>Argument</b>	'value' sets the angle, in radians.
<b>Example</b>	mySecant = sec(myAngle)
<b>Comments</b>	None.

#### 3.2.4.8 Double csc(double value)

<b>Syntax</b>	<b>Double csc(double value)</b>
<b>Description</b>	Returns the cosecant of the specified angle.
<b>Argument</b>	'value' sets the angle, in radians.
<b>Example</b>	myCosecant = csc(myAngle)
<b>Comments</b>	None.

#### 3.2.4.9 Double cot(double value)

<b>Syntax</b>	<b>Double cot(double value)</b>
<b>Description</b>	Returns the cotangent of the specified angle.
<b>Argument</b>	'value' sets the angle, in radians.
<b>Example</b>	myCotangent = cot(myAngle)
<b>Comments</b>	None.

#### 3.2.4.10 Double asin(double value)

<b>Syntax</b>	<b>Double asin(double value)</b>
<b>Description</b>	Returns the arcsine of the specified value, in radians.
<b>Argument</b>	'value' sets the sin value.
<b>Example</b>	myArcsine = asin(myAngle)
<b>Comments</b>	None.

#### 3.2.4.11 Double acos(double value)

<b>Syntax</b>	<b>Double acos(double value)</b>
<b>Description</b>	Returns the arccosine of the specified value, in radians.

<b>Argument</b>	'value' sets the cos value.
<b>Example</b>	myArcosine = acos(myAngle)
<b>Comments</b>	None.

#### 3.2.4.12 Double atan(double value)

<b>Syntax</b>	<b>Double atan(double value)</b>
<b>Description</b>	Returns the arctangent of the specified value, in radians.
<b>Argument</b>	'value' sets the tan value.
<b>Example</b>	myArctangent = atan(myAngle)
<b>Comments</b>	None.

#### 3.2.4.13 Double sinh(double value)

<b>Syntax</b>	<b>Double sinh(double value)</b>
<b>Description</b>	Returns the hyperbolic sine of the specified angle.
<b>Argument</b>	'value' sets the angle, in radians.
<b>Example</b>	myHyperSine= sinh(myAngle)
<b>Comments</b>	None.

#### 3.2.4.14 Double cosh(double value)

<b>Syntax</b>	<b>Double cosh(double value)</b>
<b>Description</b>	Returns the hyperbolic cosine of the specified angle.
<b>Argument</b>	'value' sets the angle, in radians.
<b>Example</b>	myHyperCosine= cosh(myAngle)
<b>Comments</b>	None.

#### 3.2.4.15 Double tanh(double value)

<b>Syntax</b>	<b>Double tanh(double value)</b>
<b>Description</b>	Returns the hyperbolic tangent of the specified angle.
<b>Argument</b>	'value' sets the angle, in radians.
<b>Example</b>	myHyperTan = tanh(myAngle)
<b>Comments</b>	None.

#### 3.2.4.16 Double abs(double value)

<b>Syntax</b>	<b>Double abs(double value)</b>
<b>Description</b>	Returns the absolute value of the specified value.



<b>Argument</b>	'value' sets the numeric value.
<b>Example</b>	myAbsolute = abs(myValue)
<b>Comments</b>	None.

#### 3.2.4.17 Double sqrt(double value)

<b>Syntax</b>	<b>Double sqrt(double value)</b>
<b>Description</b>	Returns the square root of the specified value
<b>Argument</b>	'value' sets the numeric value.
<b>Example</b>	myRoot = sqrt(myValue)
<b>Comments</b>	Numeric value must be positive or 0.

#### 3.2.4.18 Int ceil(double value)

<b>Syntax</b>	<b>Int ceil(double value)</b>
<b>Description</b>	Returns the smallest integer greater than the specified value.
<b>Argument</b>	'value' sets the numeric value.
<b>Example</b>	myInteger = ceil(myValue)
<b>Comments</b>	None.

#### 3.2.4.19 Int floor(double value)

<b>Syntax</b>	<b>Int floor(double value)</b>
<b>Description</b>	Returns the largest integer smaller than the specified value.
<b>Argument</b>	'value' sets the numeric value.
<b>Example</b>	myInteger = floor(myValue)
<b>Comments</b>	None.

#### 3.2.4.20 Double exp(double value)

<b>Syntax</b>	<b>Double exp(double value)</b>
<b>Description</b>	Returns e risen to a power indicated by value.
<b>Argument</b>	'value' sets the numeric value.
<b>Example</b>	myPower = exp(myValue)
<b>Comments</b>	None.

#### 3.2.4.21 Double log10(double value)

<b>Syntax</b>	<b>Double log10(double value)</b>
<b>Description</b>	Returns the decimal logarithmic of the specified value.

<b>Argument</b>	'value' sets the numeric value for the logarithm.
<b>Example</b>	myDecLog = log10(myValue)
<b>Comments</b>	None.

#### 3.2.4.22 Double log(double value1, [double value2])

<b>Syntax</b>	<b>Double log(double value1, [double value2])</b>
<b>Description</b>	Returns the logarithmic of the specified value, in the specified base.
<b>Argument</b>	'value1' sets the numeric value for the logarithm. 'value2' sets the base of the logarithm. If omitted, base is 'e'.
<b>Example</b>	myLog = log(myValue,myBase)
<b>Comments</b>	None.

#### 3.2.4.23 Double max(double value1, double value2)

<b>Syntax</b>	<b>Double max(double value1, double value2)</b>
<b>Description</b>	Returns the largest value of the two specified values.
<b>Argument</b>	'value1' sets the first numeric value. 'value2' sets the second numeric value.
<b>Example</b>	myMaximum = max(a,b)
<b>Comments</b>	None.

#### 3.2.4.24 Double min(double value1, double value2)

<b>Syntax</b>	<b>Double min(double value1, double value2)</b>
<b>Description</b>	Returns the smallest value of the two specified values.
<b>Argument</b>	'value1' sets the first numeric value. 'value2' sets the second numeric value.
<b>Example</b>	myMinimum = min(a,b)
<b>Comments</b>	None.

#### 3.2.4.25 Double pow(double value1, double value2)

<b>Syntax</b>	<b>Double pow(double value1, double value2)</b>
<b>Description</b>	Returns the power of the value in the specified base.
<b>Argument</b>	'value1' sets the the base. 'value2' sets the power.
<b>Example</b>	myValue = pow(myBase,myPower)
<b>Comments</b>	None.

### 3.2.4.26 Double round(double value1, int value2)

<b>Syntax</b>	<b>Double round(double value1, int value2)</b>
<b>Description</b>	Returns a rounded value.
<b>Argument</b>	'value1' sets the numeric value. 'value2' sets the number of decimal digits.
<b>Example</b>	myRounded = round(myValue,2)
<b>Comments</b>	None.

## 3.2.5 Date and Time functions

### 3.2.5.1 DateTime now()

<b>Syntax</b>	<b>DateTime now()</b>
<b>Description</b>	Returns the actual date and time.
<b>Argument</b>	None.
<b>Example</b>	thisMoment = now()
<b>Comments</b>	The default format is described as follow: H: hour (24), h: hour (12), m: minutes, s: seconds, t: AM/PM d: day, M: month, y: year. 'MM/dd/yyyy hh:mm:ss tt'

### 3.2.5.2 DateTime today()

<b>Syntax</b>	<b>DateTime today()</b>
<b>Description</b>	Returns the actual date.
<b>Argument</b>	None.
<b>Example</b>	thisDay = today()
<b>Comments</b>	None.

### 3.2.5.3 DateTime mindate()

<b>Syntax</b>	<b>DateTime mindate()</b>
<b>Description</b>	Returns the minimum value for the date (00:00:00.0000000, 1 <sup>st</sup> of January, 0001).
<b>Argument</b>	None.
<b>Example</b>	minMoment = mindate()
<b>Comments</b>	None.

### 3.2.5.4 *DateTime maxdate()*

<b>Syntax</b>	<b>DateTime maxdate()</b>
<b>Description</b>	Returns the maximum value for the date (23:59:59.9999999, 31 <sup>st</sup> of December, 9999).
<b>Argument</b>	None.
<b>Example</b>	maxMoment = maxdate()
<b>Comments</b>	None.

### 3.2.5.5 *String monthname()*

<b>Syntax</b>	<b>String monthname()</b>
<b>Description</b>	Returns the name of the actual month.
<b>Argument</b>	None.
<b>Example</b>	thisMonth = monthname()
<b>Comments</b>	None.

### 3.2.5.6 *DateTime adddays(dateTime date, int value)*

<b>Syntax</b>	<b>DateTime adddays(dateTime <i>date</i>, int <i>value</i>)</b>
<b>Description</b>	Adds the specified amount of days to the specified date.
<b>Argument</b>	'date' represents the original date. 'value' sets the number of days to be added, can be negative
<b>Example</b>	myDate=today() newDate=adddays(myDate,15)
<b>Comments</b>	None.

### 3.2.5.7 *DateTime addmonths(dateTime date, int value)*

<b>Syntax</b>	<b>DateTime addmonths(dateTime <i>date</i>, int <i>value</i>)</b>
<b>Description</b>	Adds the specified amount of months to the specified date
<b>Argument</b>	'date' represents the original date. 'value' sets the number of months to be added, can be negative
<b>Example</b>	myDate=today() newDate=addmonths(myDate,6)
<b>Comments</b>	None.

### 3.2.5.8 *DateTime addyears(dateTime date, int value)*

<b>Syntax</b>	<b>DateTime addyears(dateTime <i>date</i>, int <i>value</i>)</b>
<b>Description</b>	Adds the specified amount of years to the specified date.

<b>Argument</b>	'date' represents the original date. 'value' sets the number of years to be added, it can be negative.
<b>Example</b>	addyears(myTime,1)
<b>Comments</b>	None.

### 3.2.5.9 DateTime addhours(dateTime time, int value)

<b>Syntax</b>	<b>DateTime addhours(dateTime time, int value)</b>
<b>Description</b>	Adds the specified amount of hours to the specified time
<b>Argument</b>	'time' represents the original time. 'value' sets the number of hours to be added, it can be negative.
<b>Example</b>	addhours(myTime,5)
<b>Comments</b>	None.

### 3.2.5.10 DateTime addminutes(dateTime time, int value)

<b>Syntax</b>	<b>DateTime addminutes(dateTime time, int value)</b>
<b>Description</b>	Adds the specified amount of minutes to the specified time
<b>Argument</b>	'time' represents the original time. 'value' sets the number of minutes to be added it can be negative.
<b>Example</b>	addminutess(myTime,20)
<b>Comments</b>	None.

### 3.2.5.11 DateTime Addseconds(dateTime time, int value)

<b>Syntax</b>	<b>DateTime Addseconds(dateTime time, int value)</b>
<b>Description</b>	Adds the specified amount of seconds to the specified time.
<b>Argument</b>	'time' represents the original time. 'value' sets the number of seconds to be added, it can be negative.
<b>Example</b>	addseconds(myTime,200)
<b>Comments</b>	None.

## 3.2.6 String functions

### 3.2.6.1 *String concat(string value1, string value2, string value3...)*

<b>Syntax</b>	<b>String concat(string value1, string value2, string value3...)</b>
<b>Description</b>	Concatenates the given strings
<b>Argument</b>	Value# represents the strings to be added.
<b>Example</b>	myString = concat(string1,string2,string3)
<b>Comments</b>	The same effect is obtained by using the concatenate strings binary operator &. Example: myString = string1 & string2 & string3

### 3.2.6.2 *String formatnum(double value, string format)*

<b>Syntax</b>	<b>String formatnum(double value, string format)</b>
<b>Description</b>	Returns a string with the value specified in the format specified.
<b>Argument</b>	'value' represents the number to be formatted. 'format' sets the format of the number. For example "#.000".
<b>Example</b>	myString = formatnum(myValue, "#.000")
<b>Comments</b>	If value is 0.00123456 '#.0000': 0.0012 'e3': 1.235e-3

### 3.2.6.3 *String formatdate(dateTime value, string format)*

<b>Syntax</b>	<b>String formatdate(dateTime value, string format)</b>
<b>Description</b>	Returns the datetime specified in the format specified.
<b>Argument</b>	'value' represents the datetime to be formatted. 'format' sets the format of the date-time. For example "hh:mm:ss".
<b>Example</b>	myString = formatdate(myDateTime, "hh:mm:ss")
<b>Comments</b>	The format is described as follows: H: hour (24), h: hour (12), m: minutes, s: seconds, t: AM/PM d: day, M: month, y: year Example: 'hh:mm:ss tt dddd/dd/MM/yyyy'

## 3.2.7 Explicit conversion functions

### 3.2.7.1 *Double todouble(object value)*

<b>Syntax</b>	<b>Double todouble(object <i>value</i>)</b>
<b>Description</b>	Converts the specified value to a double.
<b>Argument</b>	'value' represents the object to be converted.
<b>Example</b>	myDouble = todouble("0.12345")
<b>Comments</b>	None.

### 3.2.7.2 *Int toint(object value)*

<b>Syntax</b>	<b>Int toint(object <i>value</i>)</b>
<b>Description</b>	Converts the specified value to an integer.
<b>Argument</b>	'value' represents the object to be converted.
<b>Example</b>	myInteger = toint(myDouble) myInteger = toint("230")
<b>Comments</b>	None.

### 3.2.7.3 *Long tolong(object value)*

<b>Syntax</b>	<b>Long tolong(object <i>value</i>)</b>
<b>Description</b>	Converts the specified value to a long integer,
<b>Argument</b>	'value' represents the object to be converted.
<b>Example</b>	myLong = tolong(myInteger) myLong = tolong("230")
<b>Comments</b>	None.

### 3.2.7.4 *Unsigned int touint(object value)*

<b>Syntax</b>	<b>Unsigned int touint(object <i>value</i>)</b>
<b>Description</b>	Converts the specified value to an unsigned integer.
<b>Argument</b>	'value' represents the object to be converted.
<b>Example</b>	myUInt = touint(myInteger)
<b>Comments</b>	None.

### 3.2.7.5 *Unsigned long toulong(object value)*

<b>Syntax</b>	<b>Unsigned long toulong(object <i>value</i>)</b>
<b>Description</b>	Converts the specified value to an unsigned long integer.

<b>Argument</b>	'value' represents the object to be converted.
<b>Example</b>	myULong = toulong(myInteger)
<b>Comments</b>	None.

### 3.2.7.6 *DateTime todatetime(object value)*

<b>Syntax</b>	<b>DateTime todatetime(object <i>value</i>)</b>
<b>Description</b>	Converts the specified value to a date Time
<b>Argument</b>	'value' represents the object to be converted.
<b>Example</b>	myDateTime = todatetime("10/10/2004 08:00:00 AM")
<b>Comments</b>	None.

### 3.2.7.7 *String tostring(object value)*

<b>Syntax</b>	<b>String tostring(object <i>value</i>)</b>
<b>Description</b>	Converts the specified value to a string.
<b>Argument</b>	'value' represents the object to be converted.
<b>Example</b>	myString = tostring(myInteger)
<b>Comments</b>	None.



# 4 BOSA SCPI Commands

The commands used by the BOSA are defined according to the Standard Commands for Programmable Instruments (SCPI) and are arranged in a command tree, also called subsystems. Every subsystem contains all commands belonging to a specific topic.

Some conventions related to the SCPI commands are:

- The part of the command shown in uppercase represents the short form of the command. The commands are case insensitive.
- Values to be input are indicated by angle brackets (<>) and are separated from the command by a colon as shown in the command syntax
- Optional values and portions of syntax are indicated by square brackets ([ ]).
- The bar (|) shows an either–or choice of data, for example a|b means either a or b, but not both simultaneously.

## 4.1 BOSA SCPI Command Summary

The following table gives an overview of the command tree. The commands and a page reference for more detailed information are also shown. Some commands are application specific and only work if the application is running. Their availability depends on your BOSA option (See section 2 in the User's Guide). The following abbreviations are used to indicate the command applicability: BOSA (B), BOSA PHASE (BP), TUNABLE LASER (TLS), COMPONENT ANALYZER (CA).

INSTRUMENT SUBSYSTEM COMMANDS			
INSTRument:STATe:MODE ?	Queries the current application.	B,TLS,CA	42
INSTRument:STATe:MODE <MAIN BOSA TLS CA>	Sets the application. To switch between applications you must always start from the MAIN	B,TLS,CA	42

	window.		
INSTrument:STATe:<HOLD? RUN?>	Queries the current measurement state, and reads whether the BOSA is on.	B,CA	42
INSTrument:STATe:<HOLD RUN> <1 0>	Sets the BOSA's state: HOLD or RUN.	B, CA	42
<b>DISPLAY SUBSYSTEM COMMANDS</b>			
DISPlay[:WINDow]:TRACe:Y[:SCALe]:AUTO [ONCE]	Sets an autoscale in Y axis.	B,CA	43
DISPlay[:WINDow]:TRACe:Y[:SCALe]:BOTTOm <value> [DBM MW]	Sets the lowest value of the vertical axis.	B,CA	43
DISPlay[:WINDow]:TRACe:Y[:SCALe]:BOTTOm?	Queries the lowest value for of the vertical axis.	B,CA	43
DISPlay[:WINDow]:TRACe:Y[:SCALe]:PDIVision <numeric_value> [DBM MW]	Sets power per division of the vertical axis.	B,CA	43
DISPlay[:WINDow]:TRACe:Y[:SCALe]:PDIVision?	Queries power per division for the vertical axis.	B,CA	43
DISPlay[:WINDow]:TRACe:Y[:SCALe]:RLEVel <numeric_value> [DBM MW]	Sets the reference level value in the vertical axis.	B,CA	44
DISPlay[:WINDow]:TRACe:Y[:SCALe]:RLEVel?	Queries the reference level value in the vertical axis.	B,CA	44
DISPlay[:WINDow]:TRACe:Y[:SCALe]:NORMAlize <0 1 OFF ON>	Turns the normalization of the Y scale on or off	CA	44
DISPlay[:WINDow]:TRACe:Y[:SCALe]: NORMAlize?	Queries the normalization of the Y scale.	CA	44
DISPlay[:WINDow]:TRACe:Y:SPACing <LOGarithmic LINEar>	Sets the vertical scale to logarithmic or lineal.	B	44
DISPlay[:WINDow]:TRACe:Y:SPACing?	Queries the current vertical scale.	B	45
DISPlay[:WINDow]:TRACe:X WAVElength FREQuency	Changes the horizontal axis to the specified units.	B,CA	45
DISPlay[:WINDow]:TRACe:X?	Returns the current units of the horizontal axis.	B,CA	45
DISPlay[:WINDow]:TRACe:STATe <A B M1 M2 M3 M4>,<ON OFF>	Change the active trace	B,CA	45

DISPlay[:WINDow]:TRACe:STATe?	Queries de active trace	B, CA	45
DISPlay[:WINDow]:GRAPHICSElection <C+L   O>	Changes the active graphic.	B	45
DISPlay[:WINDow]:GRAPHICSElection ?>	Queries the active graphic.	B	46
DISPlay[:WINDow]:GRAPHICVIEW <C+L   O   O+C+L>	Expands to full screen the desired graphic.	B	46
DISPlay[:WINDow]:GRAPHICVIEW ?>	Queries the displayed graphic.	B	46
<b>SENSE SUBSYSTEM COMMANDS</b>			
SENSe:WAVelength:CENTer <value> [NM PM GHZ THZ]	Sets the center value of the horizontal axis.	B,CA	46
SENSe:WAVelength:CENTer?	Queries the center value of the horizontal axis.	B,CA	47
SENSe:WAVelength:SINGLE <1 0 ON OFF>	Set a single or continuous sweep	TLS	47
SENSe:WAVelength:SINGLE?	Queries the type of sweeping	TLS	47
SENSe:WAVelength:SMOOTH <numeric_value> [NM PM GHZ THZ]	Performs a smooth operation over the active trace	CA	47
SENSe:WAVelength:SMOOTH?	Queries the current smoothing value	CA	47
SENSe:WAVelength:SPAN <numeric_value> [NM PM GHZ THZ]	Sets the span value of the horizontal axis.	B,CA	48
SENSe:WAVelength:SPAN?	Queries the span value of the horizontal axis.	B,CA	48
SENSe:WAVelength:SPEED <numeric_value> [NM PM GHZ THZ]	Sets the sweep speed for the internal tunable laser	TLS	48
SENSe:WAVelength:SPEED?	Queries the current sweep speed of the internal tunable laser.	TLS	48
SENSe:WAVelength:STATic <numeric_value> [NM PM GHZ THZ]	Sets the output wavelength for the internal tunable laser	TLS	48
SENSe:WAVelength:STATic?	Queries the current	TLS	48

	tuned wavelength for the internal tunable laser.		
SENSe:WAVelength:StARt <numeric_value> [NM PM GHZ THZ]	Sets the minimum (left) value of the horizontal axis or the starting wavelength for a sweep.	B,CA,TLS	49
SENSe:WAVelength:StARt?	Queries the minimum (left) value of the horizontal axis or the starting wavelength for a sweep.	B,CA,TLS	49
SENSe:WAVelength:StOP <numeric_value> [NM PM GHZ THZ]	Sets the maximum (right) value of the horizontal axis or the stopping wavelength for a sweep.	B,CA,TLS	49
SENSe:WAVelength:StOP?	Queries the maximum (right) value of the horizontal axis or the stopping wavelength for a sweep.	B,CA,TLS	49
SENSe:WAVelength:RESolution <numeric value>	Sets the BOSA resolution bandwidth	B	49
SENSe:WAVelength:RESolution?	Queries the BOSA sampling rate.	B	50
SENSe:WAVelength:SMODE <HR HS>	Sets the measurement speed mode	CA	50
SENSe:WAVelength:SMODE?	Queries the measurement speed mode	CA	50
SENSe:AVERage:COUnT <4 8 12 32 CONT>	Sets the number of averages.	B,CA	50
SENSe:AVERage:COUnT?	Queries the current count of averages.	B,CA	50
SENSe:AVERage:STATe <0 1 OFF ON>	Turns the averaging mode on or off.	B,CA	51
SENSe:AVERage:STATe?	Queries the state of averaging mode.	B,CA	51
SENSe:AVERage:CORRelate <0 1 ON OFF>	Enables/disables the correlation function.	B,CA	51
SENSe:AVERage:CORRelate?	Queries the state of	B,CA	51

	the correlation function.		
SENSe:AVERage:CORRelate:CENTer <numeric_value> [NM]	Sets the center of the correlation zone on the active trace.	B,CA	51
SENSe:AVERage:CORRelate:CENTer?	Queries the center of the correlation zone on the active trace.	B,CA	52
SENSe:AVERage:CORRelate:SPAN <numeric_value> [NM]	Sets the width of the correlation zone on the active trace.	B,CA	52
SENSe:AVERage:CORRelate:SPAN?	Queries the width of the correlation zone on the active trace.	B,CA	52
SENSe:NOISezeroing	Performs a noise zeroing	B,CA	52
SENSe:SWITCH <0 1 ON OFF>	Switches on or off the internal tunable laser	TLS	52
SENSe:SWITCH?	Queries if the internal tunable laser is switched on or off	TLS	53
SENSe:SWEEP <0 1 ON OFF>	Starts or stops a sweep	TLS	53
SENSe:SWEEP?	Queries the sweeping status	TLS	53
SENSe:LASer <C+L   O>	Selects the internal laser to operate with through the TLS OUT front panel port.	TLS	53
SENSe:LASer?	Queries the actual internal laser that is being externally used for user purposes.	TLS	53
SENSe:BAND <C+L   O   O+C+L>	Change the measurement band.	B,CA	54
SENSe:BAND?>	Queries the actual measured band.	B,CA	54
<b>INPUT SUBSYSTEM COMMANDS</b>			
INPut:SPARameters <IL RL IL&RL>	Sets the type of measurement to display	CA	54
INPut:SPARameters?	Queries the type of measurement being displayed	CA	54
INPut:POLarization <1+2 1 2 1&2>	Sets the polarization	B	54

	option measurement		
INPut:POLarization <PDL MAX MIN SIMUL>	Sets the displayed measure	CA (with optX30)	55
INPut:POLarization <INDEP 1 2 SIMUL>	Sets the displayed measure	CA (without optX30)	55
INPut:POLarization?	Queries the measured polarization	B,CA	55
INPut:POLarization:MUELLermode <ON OFF 1 0>	Sets the Mueller Mode on or off	CA (with optX30)	55
INPut:POLarization:MUELLermode?	Queries if the Mueller Mode is active	CA (with optX30)	56
INPut:POWer?	Queries the total SUT power	B	56
<b>CALCULATE SUBSYSTEM COMMANDS</b>			
CALCulate:MARKer:AOff	Disable the marker on the active trace.	B,CA	56
CALCulate:MARKer:STATe <ON OFF 1 0>	Enables/disables marker on the active trace.	B,CA	56
CALCulate:MARKer:STATe?	Queries the state of the marker on the active trace.	B,CA	57
CALCulate:MARKer:MODE <TRCK FIXX FIXXY>	Sets the marker' behaviour.	B, CA	57
CALCulate:MARKer:MODE?	Queries the marker's behaviour.	B, CA	57
CALCulate:MARKer:MAXimum	Moves the marker on the active trace to the maximum power level of the trace.	B,CA	57
CALCulate:MARKer:MAXimum:NEXT	Moves the marker to the next power maximum.	B,CA	57
CALCulate:MARKer:MAXimum:RIGHT	Moves the marker to the next power maximum to the right.	B,CA	58
CALCulate:MARKer:MAXimum:LEFT	Moves the marker to the next power maximum to the left.	B,CA	58

CALCulate:MARKer:SCEN	Moves the marker on the active trace to the center of the screen.	B,CA	58
CALCulate:MARKer:X <numeric_value> [NM PM GHZ THZ]	Moves the marker on the active trace to the specified point on the horizontal axis.	B,CA	58
CALCulate:MARKer:X?	Queries the horizontal coordinate of the marker on the active trace.	B,CA	58
CALCulate:MARKer:Y <numeric_value> [DBM MW]	Moves the marker on the active trace to the specified point on the vertical axis.	B,CA	59
CALCulate:MARKer:Y?	Queries the vertical coordinate of the marker on the active trace.	59	59
CALCulate:MARKer:THREshold <numeric_value> [DB]	Sets the peak excursion value for the peak search functions.	59	59
CALCulate:MARKer:THREshold?	Queries the current peak excursion value for the peak search functions.	B,CA	59
CALCulate:MARKer:READout FREQuency   WAVElengt	Sets the horizontal units of the marker	B,CA	59
CALCulate:MARKer:READout?	Queries the horizontal units of the marker	B,CA	60
CALCulate:MARKer:SRLevel	Moves the reference level to the vertical value of the marker on the active trace.	B,CA	60
CALCulate:MARKer:POL?	In polarimetry, queries the S parameters (S1, S2, S3) of the marker	B	60
CALCulate:MARKer:FUNCtion:DELTA [:STATe] <ON OFF 0 1>	Turns the delta marker function on or off.	B,CA	60

CALCulate:MARKer:FUNCTION:DELTA[:STATE]?	Queries the state of the delta marker function	B,CA	60
CALCulate:MARKer:FUNCTION:DELTA:RESet	Sets the marker offset to the current position of the marker	B,CA	60
CALCulate:MARKer:FUNCTION:DELTA:X:OFFSet?	Queries the horizontal offset between the marker and its reference	B,CA	61
CALCulate:MARKer:FUNCTION:DELTA:X:REFerece?	Queries the horizontal value of the reference marker	B,CA	61
CALCulate:MARKer:FUNCTION:DELTA:Y:OFFSet?	Queries the vertical offset between the marker and its reference	B,CA	61
CALCulate:MARKer:FUNCTION:DELTA:Y:REFerece?	Queries the horizontal value of the reference marker	B,CA	61
CALCulate:MARKer:FUNCTION:DELTA:POL?	In polarimetry, queries the S parameters (S1, S2, S3) of the offset	B	61
CALCulate:MARKer:FUNCTION:DELTA:ANG?	In polarimetry, queries the S angle between marker and offset	B	61
CALCulate:MAXimum[:STATE] <0 1 ON OFF>	Enables/disables the Max Hold function.	B,CA	61
CALCulate:MAXimum[:STATE]?	Queries the state of the Max Hold function.	B,CA	62
CALCulate:MINimum[:STATE] ON OFF	Enables/disables the Min Hold function.	B,CA	62
CALCulate:MINimum[:STATE]?	Queries the state of the Min Hold function.	B,CA	62
CALCulate:TPOWER ON OFF	Enables/disables the power integral function	B,CA	63
CALCulate:TPOWER:[DATA]?	Returns the power integral.	B,CA	63



CALCulate:TPOWer:IRANge:UPPer<numeric_value> [NM PM GHZ THZ]	Sets the upper bound for power integration.	B,CA	63
CALCulate:TPOWer:IRANge:UPPer?	Returns the upper bound for power integration.	B,CA	63
CALCulate:TPOWer:IRANge:LOWer <numeric_value> [NM PM GHZ THZ]	Sets the lower bound for power integration.	B,CA	63
CALCulate:TPOWer:IRANge:LOWer?	Returns the lower bound for power integration.	B,CA	63
CALCulate:AUXINput:POWer?	Returns the power meter measurement from the Aux Input port.	TLS	64
<b>TRACE SUBSYSTEM COMMANDS</b>			
TRACe[:DATA]:COUNT?	Queries the number of points displayed on the screen.	B,CA	64
TRACe[:DATA]?	Queries the entire spectrum measured by the BOSA.	B,CA	64
TRACe[:DATA]:MAXimum:X?	Queries the X coordinate of the absolute maximum.	B,CA	65
TRACe[:DATA]:MAXimum:Y?	Queries the Y coordinate of the absolute maximum.	B,CA	65
FORMat[:DATA] <ASCII REAL>,<length>	Sets the output data format for queries either to ASCII with a given length or REAL (64 bits).	B,CA	65
FORMat[:DATA]?	Returns the current output data format of the GPIB queries.	B,CA	65
<b>MMEMORY SUBSYSTEM COMMANDS</b>			
MMEMorY:STORe:TRACe <name>.<type>,[ink_saving_mode]	Stores the traces displayed in a file in the specified format.	B,CA	66
MMEMorY:DELete:TRACe <name>	Deletes the specified file	B,CA	66
MMEMorY:LOAD:TRACe <M1 M2 M3 M4>,<name>	Loads the specified file	B,CA	66

## 4.2 BOSA SCPI Command Definitions

### 4.2.1 INSTRument Subsystem Commands

#### 4.2.1.1 *INSTRument:STATe:MODE ?*

<b>Syntax</b>	<b>INSTRument:STATe:MODE ?</b>
<b>Description</b>	Queries the current application.
<b>Response</b>	'MAIN' if the unit is running the main menu, and <BOSA TLS CA> if it is running one of the applications.
<b>Example</b>	INST:STAT:MODE ?

#### 4.2.1.2 *INSTRument:STATe:MODE <MAIN|BOSA|TLS|CA>*

<b>Syntax</b>	<b>INSTRument:STATe:MODE &lt;MAIN BOSA TLS CA&gt;</b>
<b>Description</b>	Sets the application. To switch between applications you must always start from the MAIN window.
<b>Example</b>	INST:STAT:MODE BOSA

#### 4.2.1.3 *INSTRument:STATe:<HOLD?|RUN?>*

<b>Syntax</b>	<b>INSTRument:STATe:&lt;HOLD? RUN?&gt;</b>
<b>Description</b>	Queries the current measurement state, and reads whether the BOSA is on.
<b>Response</b>	'ON' if it is running, and 'OFF' if it is stopped.
<b>Example</b>	INST:STAT: RUN?

#### 4.2.1.4 *INSTRument:STATe:<HOLD|RUN>*

<b>Syntax</b>	<b>INSTRument:STATe:&lt;HOLD RUN&gt; &lt;1 0&gt;</b>
<b>Description</b>	Sets the BOSA's state: HOLD or RUN.
<b>Argument</b>	'HOLD' to freeze the trace on the screen, the BOSA continues working on a second plane and 'RUN' to start and stop the application.
<b>Example</b>	INST:STAT:RUN 1

## 4.2.2 DISPlay Subsystem Commands

### 4.2.2.1 *DISPlay[:WINDow]:TRACe:Y[:SCALe]:AUTO [ONCE]*

Syntax	DISPlay[:WINDow]:TRACe:Y[:SCALe]:AUTO [ONCE]
Description	Performs an auto scale in the vertical axis.
Example	DISP:TRAC:Y:AUTO

### 4.2.2.2 *DISPlay[:WINDow]:TRACe:Y[:SCALe]:BOTTom*

Syntax	DISPlay[:WINDow]:TRACe:Y[:SCALe]:BOTTom <value> [DBM MW]
Description	Sets the lowest value of the vertical axis.
Argument	'value' is the lowest value of the vertical axis. It can be a numeric value between 20 and -80 dBm or the Ref Level, using the key 'rlev'. Units can be 'DBM' for logarithmic scale or 'MW' for linear scale.
Example	DISP:TRAC:Y:BOTT -70 DBM

### 4.2.2.3 *DISPlay[:WINDow]:TRACe:Y[:SCALe]:BOTTom?*

Syntax	DISPlay[:WINDow]:TRACe:Y[:SCALe]:BOTTom?
Description	Queries the lowest value of the vertical axis.
Response	The lowest value of the vertical axis, expressed in the current units.
Example	DISP:TRAC:Y:BOTT?

### 4.2.2.4 *DISPlay[:WINDow]:TRACe:Y[:SCALe]:PDIVision*

Syntax	DISPlay[:WINDow]:TRACe:Y[:SCALe]:PDIVision <numeric_value> [DBM MW]
Description	Sets power per division of the vertical axis.
Argument	'numeric_value' is the power per division (max. 10 dB) for the vertical axis. Units are optional. Values must match the Y scale units.
Example	DISP:TRAC:Y:PDIV 10 DBM

### 4.2.2.5 *DISPlay[:WINDow]:TRACe:Y[:SCALe]:PDIVision?*

Syntax	DISPlay[:WINDow]:TRACe:Y[:SCALe]:PDIVision?
Description	Queries power per division for the vertical axis.

<b>Response</b>	The power per division for the vertical axis, expressed in the current units.
<b>Example</b>	DISP:TRAC:Y:PDIV?

#### 4.2.2.6 *DISPlay[:WINDow]:TRACe:Y[:SCALe]:RLEVel*

<b>Syntax</b>	<b>DISPlay[:WINDow]:TRACe:Y[:SCALe]:RLEVel &lt;numeric_value&gt; [DBM MW]</b>
<b>Description</b>	Sets the reference level value in the vertical axis.
<b>Argument</b>	'numeric_value' is the value for the reference level.
<b>Example</b>	DISP:TRAC:Y:RLEV 0 DBM

#### 4.2.2.7 *DISPlay[:WINDow]:TRACe:Y[:SCALe]:RLEVel?*

<b>Syntax</b>	<b>DISPlay[:WINDow]:TRACe:Y[:SCALe]:RLEVel?</b>
<b>Description</b>	Queries the reference level value in the vertical axis.
<b>Response</b>	The vertical value for the reference level, expressed in the current units.
<b>Example</b>	DISP:TRAC:Y:RLEV?

#### 4.2.2.8 *DISPlay[:WINDow]:TRACe:Y[:SCALe]:NORMalize*

<b>Syntax</b>	<b>DISPlay[:WINDow]:TRACe:Y[:SCALe]:NORMalize &lt;0 1 OFF ON&gt;</b>
<b>Description</b>	Turns the normalization of the Y scale on or off.
<b>Example</b>	DISP:TRAC:Y:NORM 0

#### 4.2.2.9 *DISPlay[:WINDow]:TRACe:Y[:SCALe]:NORMalize?*

<b>Syntax</b>	<b>DISPlay[:WINDow]:TRACe:Y[:SCALe]:NORMalize?</b>
<b>Description</b>	Queries the normalization of the Y scale.
<b>Response</b>	It returns 1 or On if the normalization is applied to the Y scale.
<b>Example</b>	DISP:TRAC:Y:NORM?

#### 4.2.2.10 *DISPlay[:WINDow]:TRACe:Y:SPACing*

<b>Syntax</b>	<b>DISPlay[:WINDow]:TRACe:Y:SPACing &lt;LOGarithmic LINear&gt;</b>
<b>Description</b>	Sets the vertical scale.
<b>Argument</b>	It can be 'logarithmic' or 'linear' for the vertical axis.
<b>Example</b>	DISP:TRAC:Y:SPAC LIN

#### 4.2.2.11 *DISPlay[:WINDow]:TRACe:Y:SPACing?*

<b>Syntax</b>	<b>DISPlay[:WINDow]:TRACe:Y:SPACing?</b>
<b>Description</b>	Queries the vertical scale.
<b>Response</b>	It returns the current spacing of the Y axis: 'Spacing Logarithmic' or 'Spacing Linear'.
<b>Example</b>	DISP:TRAC:Y:SPAC?

#### 4.2.2.12 *DISPlay[:WINDow]:TRACe:X WAVelength|FREQuency*

<b>Syntax</b>	<b>DISPlay[:WINDow]:TRACe:X WAVelength FREQuency</b>
<b>Description</b>	Changes the x axis to the specified units, nm or GHz
<b>Example</b>	DISP:TRAC:X WAV

#### 4.2.2.13 *DISPlay[:WINDow]:TRACe:X?*

<b>Syntax</b>	<b>DISPlay[:WINDow]:TRACe:X?</b>
<b>Description</b>	Returns the current units of the x axis.
<b>Example</b>	DISP:TRAC:X?

#### 4.2.2.14 *DISPlay[:WINDow]:TRACe:STATe*

<b>Syntax</b>	<b>DISPlay[:WINDow]:TRACe:STATe &lt;M1 M2 M3 M4&gt;,&lt;ON OFF&gt;</b>
<b>Description</b>	Changes the active trace.
<b>Example</b>	DISP:TRAC:STAT M2,ON
<b>Comments</b>	If trace M2 is activated and it does not exist, the Active trace is copied in trace M2. If the trace exist it becomes the active trace. The OFF state clear the selected trace

#### 4.2.2.15 *DISPlay[:WINDow]:TRACe:STATe?*

<b>Syntax</b>	<b>DISPlay[:WINDow]:TRACe:STATe?</b>
<b>Description</b>	Queries the active trace.
<b>Example</b>	DISP:TRAC:STAT?

#### 4.2.2.16 *DISPlay[:WINDow]:GRAPHICSElection <C+L|O>*

<b>Syntax</b>	<b>DISPlay[:WINDow]:GRAPHICSElection &lt;C+L O&gt;</b>
---------------	--

<b>Description</b>	Changes the active graphic to C+L or O band graphic. The sent commands afterwards are enabled over the active graphic. (Only available for BOSA O+C+L)
<b>Example</b>	DISP:GRAPHSEL C+L

#### 4.2.2.17 DISPlay[:WINDow]:GRAPHICSElection ?

<b>Syntax</b>	DISPlay[:WINDow]:GRAPHICSElection <C+L O>
<b>Description</b>	Queries the active graphic. (Option only valid for BOSA O+C+L).
<b>Response</b>	'C+L' if the selected graph corresponds to the C+L band, 'O' if the selected graph corresponds to the O band.
<b>Example</b>	DISP:GRAPHSEL ?

#### 4.2.2.18 DISPlay[:WINDow]:GRAPHICVIEW<C+L|O|O+C+L>

<b>Syntax</b>	DISPlay[:WINDow]:GRAPHICVIEW <C+L O O+C+L>
<b>Description</b>	Expands to full screen the desired graphic, when O+C+L option is sent, screen is divided into 2 graphs. O+C+L option is only available when simultaneous bands are measured, see command 4.2.3.38. (Only available for BOSA O+C+L)
<b>Example</b>	DISP:GRAPHVIEW C+L

#### 4.2.2.19 DISPlay[:WINDow]:GRAPHICVIEW ?

<b>Syntax</b>	DISPlay[:WINDow]:GRAPHICVIEW ?
<b>Description</b>	Queries the displayed graphic. (Option only valid for BOSA O+C+L)
<b>Response</b>	'C+L' if the displayed graph corresponds to the C+L band, 'O' if the displayed graph corresponds to the O band, 'O+C+L' if simultaneous mode.
<b>Example</b>	DISP:GRAPHVIEW ?

### 4.2.3 SENSE Subsystem Commands

#### 4.2.3.1 SENSE:WAVelength:CENTer

<b>Syntax</b>	SENSe:WAVelength:CENTer <value> [NM PM GHZ THZ]
<b>Description</b>	Sets the center value of the horizontal axis.
<b>Argument</b>	It can be a numeric value or 'MAX'. If 'MAX' is used, the

	maximum power level displayed in the screen is moved to the center.
<b>Example</b>	SENS:WAV:CENT 1550 NM

#### 4.2.3.2 SENSE:WAVelength:CENTer?

<b>Syntax</b>	SENSe:WAVelength:CENTer?
<b>Description</b>	Queries the center value of the horizontal axis.
<b>Response</b>	The center value of the horizontal axis, expressed in nm or GHz
<b>Example</b>	SENS:WAV:CENT?

#### 4.2.3.3 SENSE:WAVelength:SINGLE

<b>Syntax</b>	SENSe:WAVelength:SINGLE 1 0 ON OFF
<b>Description</b>	Sets the type of sweeping for the internal tunable laser
<b>Argument</b>	'1' or ON set a single sweep, '0' or OFF set a continuous sweep.
<b>Example</b>	SENS:WAV:SINGLE ON

#### 4.2.3.4 SENSE:WAVelength:SINGLE?

<b>Syntax</b>	SENSe:WAVelength:SINGLE?
<b>Description</b>	Queries the current type of sweeping of the internal tunable laser
<b>Response</b>	ON means a single sweep is set, OFF means a continuous sweep is set.
<b>Example</b>	SENS:WAV:SINGLE?

#### 4.2.3.5 SENSE:WAVelength:SMOOTH

<b>Syntax</b>	SENSe:WAVelength:SMOOTH <value> [NM PM GHZ THZ]
<b>Description</b>	Performs a smoothing operation over the active trace
<b>Argument</b>	The spectral width averaged in the smoothing
<b>Example</b>	SENS:WAV:SMOOTH 1 GHZ

#### 4.2.3.6 SENSE:WAVelength:SMOOTH?

<b>Syntax</b>	SENSe:WAVelength:SMOOTH?
<b>Description</b>	Queries the current value of the smoothing operation
<b>Response</b>	The spectral width averaged in the smoothing
<b>Example</b>	SENS:WAV:SMOOTH?

#### 4.2.3.7 *SENSe:WAVelength:SPAN*

<b>Syntax</b>	<b>SENSe:WAVelength:SPAN &lt;numeric_value&gt; [NM PM GHZ THZ]</b>
<b>Description</b>	Sets the span value of the horizontal axis.
<b>Argument</b>	The span value of the horizontal axis.
<b>Example</b>	SENS:WAV:SPAN 2 NM

#### 4.2.3.8 *SENSe:WAVelength:SPAN?*

<b>Syntax</b>	<b>SENSe:WAVelength:SPAN?</b>
<b>Description</b>	Queries the span value of the horizontal axis.
<b>Response</b>	The span value of the horizontal axis, expressed in nm or GHz.
<b>Example</b>	SENS:WAV:SPAN?

#### 4.2.3.9 *SENSe:WAVelength:SPEED*

<b>Syntax</b>	<b>SENSe:WAVelength:SPEED &lt;numeric_value&gt; [NM PM GHZ THZ]</b>
<b>Description</b>	Sets the sweep speed of the internal tunable laser
<b>Argument</b>	The sweep speed. Units are given per second.
<b>Example</b>	SENS:WAV:SPEED 2 NM

#### 4.2.3.10 *SENSe:WAVelength:SPEED?*

<b>Syntax</b>	<b>SENSe:WAVelength:SPEED?</b>
<b>Description</b>	Queries the current sweep speed of the internal tunable laser.
<b>Argument</b>	The sweep speed. Units are given per second.
<b>Example</b>	SENS:WAV:SPEED?

#### 4.2.3.11 *SENSe:WAVelength:STAT*

<b>Syntax</b>	<b>SENSe:WAVelength:STAT &lt;numeric_value&gt; [NM PM GHZ THZ]</b>
<b>Description</b>	Sets the tuning wavelength for the internal tunable laser.
<b>Argument</b>	The tunable laser output wavelength
<b>Example</b>	SENS:WAV:STAT 1549 NM

#### 4.2.3.12 *SENSe:WAVelength:STAT?*

<b>Syntax</b>	<b>SENSe:WAVelength:STAT?</b>
<b>Description</b>	Queries the tuned wavelength for the internal tunable laser.
<b>Response</b>	The tuned wavelength



Example	SENS:WAV:STAT?
---------	----------------

#### 4.2.3.13 SENSE:WAVelength:STARt

Syntax	SENSe:WAVelength:STARt <numeric_value> [NM PM GHZ THZ]
Description	Sets the minimum (left) value of the horizontal axis or the starting wavelength for a sweep.
Argument	The value for the minimum in the horizontal axis or the starting wavelength for a sweep.
Example	SENS:WAV:START 1549 NM

#### 4.2.3.14 SENSE:WAVelength:STARt?

Syntax	SENSe:WAVelength:STARt?
Description	Queries the minimum (left) value of the horizontal axis or the starting wavelength for a sweep.
Response	The value for the minimum in the horizontal axis or the starting wavelength for a sweep, expressed in nanometers.
Example	SENS:WAV:START?

#### 4.2.3.15 SENSE:WAVelength:STOP

Syntax	SENSe:WAVelength:STOP <numeric_value> [NM PM GHZ THZ]
Description	Sets the maximum (right) value of the horizontal axis or the stopping wavelength for a sweep.
Argument	The value for the maximum in the horizontal axis.
Example	SENS:WAV:STOP? 1551 NM

#### 4.2.3.16 SENSE:WAVelength:STOP?

Syntax	SENSe:WAVelength:STOP?
Description	Queries the maximum (right) value of the horizontal axis or the stopping wavelength for a sweep.
Response	The value for the maximum of the horizontal axis or the stopping wavelength for a sweep, expressed in nm or GHz.
Example	SENS:WAV:STOP?

#### 4.2.3.17 SENSE:WAVelength:RESolution

Syntax	SENSe:WAVelength:RESolution <value> [NM PM GHZ THZ]
--------	---

<b>Description</b>	Sets the resolution bandwidth.
<b>Argument</b>	The spectral filter width.
<b>Example</b>	SENS:WAV:RES 1 GHZ

#### 4.2.3.18 SENSE:WAVelength:RESolution?

<b>Syntax</b>	SENSe:WAVelength:RESolution?
<b>Description</b>	Queries the resolution bandwidth.
<b>Response</b>	The spectral filter width.
<b>Example</b>	SENS:WAV:RES?

#### 4.2.3.19 SENSE:WAVelength:SMODE

<b>Syntax</b>	SENSe:WAVelength:SMODE <HR HS>
<b>Description</b>	Change the speed mode of the component analyser measurement.
<b>Argument</b>	The value 'HR' stands for the High Resolution mode while the 'HS' sets the mode to the High Speed mode. Refer to the User Guide for further explaining.
<b>Example</b>	SENS:WAV:SMOD HR

#### 4.2.3.20 SENSE:WAVelength:SMODE?

<b>Syntax</b>	SENSe:WAVelength:SMODE?
<b>Description</b>	Queries the actual measurement speed mode
<b>Response</b>	'HR' stands for the High Resolution mode while the 'HS' stands for the High Speed mode
<b>Example</b>	SENS:WAV:SMOD?

#### 4.2.3.21 SENSE:AVERage:COUNt

<b>Syntax</b>	SENSe:AVERage:COUNt 4 8 12 32 CONT
<b>Description</b>	Sets the number of averages.
<b>Argument</b>	Available number of averages: 4, 8, 12, 32 or CONTINUOUS.
<b>Example</b>	SENS:AVER:COUN CONT

#### 4.2.3.22 SENSE:AVERage:COUNt?

<b>Syntax</b>	SENSe:AVERage:COUNt?
<b>Description</b>	Queries the current number of averages.

<b>Response</b>	The current count of averages. It is recommended to use <i>wait</i> command to let the BOSA reach the number of averages set.
<b>Example</b>	SENS:AVER:COUN?

#### 4.2.3.23 SENSE:AVERage:STATe

<b>Syntax</b>	SENSe:AVERage:STATe 0 1 OFF ON
<b>Description</b>	Enables or disables the averaging mode.
<b>Argument</b>	'1' or 'ON' to start a new average, and '0' or 'OFF' to stop it.
<b>Example</b>	SENS:AVER:STAT ON

#### 4.2.3.24 SENSE:AVERage:STATe?

<b>Syntax</b>	SENSe:AVERage:STATe?
<b>Description</b>	Queries the state of averaging mode.
<b>Response</b>	The current state of the averaging mode: 'ON' if averaging, 'OFF' if stopped.
<b>Example</b>	SENS:AVER:STAT?

#### 4.2.3.25 SENSE:AVERage:CORRelate

<b>Syntax</b>	SENSe:AVERage:CORRelate 0 1 ON OFF
<b>Description</b>	Enables or disables the correlation function in order to lock the trace relatively to a peak.
<b>Argument</b>	'1' or 'on' to enable the correlation; '0' or 'off' to disable it.
<b>Example</b>	SENS:AVER:CORR ON

#### 4.2.3.26 SENSE:AVERage:CORRelate?

<b>Syntax</b>	SENSe:AVERage:CORRelate?
<b>Description</b>	Queries the state of the correlation function.
<b>Response</b>	'ON' if it is enabled, 'OFF' if not.
<b>Example</b>	SENS:AVER:CORR?

#### 4.2.3.27 SENSE:AVERage:CORRelate:CENTer

<b>Syntax</b>	SENSe:AVERage:CORRelate:CENTer <numeric_value> [NM]
<b>Description</b>	Sets the center of the correlation zone.
<b>Argument</b>	<numeric_value> is the value for the center of the correlation zone.

<b>Example</b>	SENS:AVER:CORR:CENT 1550 NM
<b>Comments</b>	Referred to the active trace.

#### 4.2.3.28 SENSE:AVERage:CORRelate:CENTer?

<b>Syntax</b>	SENSe:AVERage:CORRelate:CENTer?
<b>Description</b>	Queries the center of the correlation zone.
<b>Response</b>	The center of the correlation zone, expressed in nm or GHz.
<b>Example</b>	SENS:AVER:CORR:CENT?
<b>Comments</b>	Referred to the active trace.

#### 4.2.3.29 SENSE:AVERage:CORRelate:SPAN

<b>Syntax</b>	SENSe:AVERage:CORRelate:SPAN <numeric_value> [NM]
<b>Description</b>	Sets the width of the correlation zone.
<b>Argument</b>	<numeric_value> is the value for the span of the correlation zone.
<b>Example</b>	SENS:AVER:CORR:SPAN 0.01 NM
<b>Comments</b>	Referred to the active trace.

#### 4.2.3.30 SENSE:AVERage:CORRelate:SPAN?

<b>Syntax</b>	SENSe:AVERage:CORRelate:SPAN?
<b>Description</b>	Queries the span of the correlation zone.
<b>Response</b>	The span of the correlation zone, expressed in nm or GHz.
<b>Example</b>	SENS:AVER:CORR:SPAN?
<b>Comments</b>	Referred to the active trace.

#### 4.2.3.31 SENSE:NOISezeroing

<b>Syntax</b>	SENSe:NOISezeroing
<b>Description</b>	Performs a noise zeroing.
<b>Example</b>	SENS:NOIS
<b>Comments</b>	This function carries out an immediate improvement in the noise pattern of the current spectrum.

#### 4.2.3.32 SENSE:SWITCH

<b>Syntax</b>	SENSe:SWITCH 1 0 ON OFF
---------------	-------------------------

<b>Description</b>	Switches on or off the internal tunable laser
<b>Arguments</b>	'1' or ON switch on the laser, '0' or OFF switch off the laser.
<b>Example</b>	SENS:SWITCH ON

#### 4.2.3.33 SENSE:SWITCH?

<b>Syntax</b>	SENSe:SWITCH?
<b>Description</b>	Queries the internal tunable laser is switched on or off
<b>Response</b>	ON if the laser switched on, OFF if the laser is switched off.
<b>Example</b>	SENS:SWITCH?

#### 4.2.3.34 SENSE:SWEEP

<b>Syntax</b>	SENSe:SWEEP 1 0 ON OFF
<b>Description</b>	Start or stops a sweep
<b>Arguments</b>	'1' or ON starts the sweep, '0' or OFF stops the sweep.
<b>Example</b>	SENS:SWEEP ON

#### 4.2.3.35 SENSE:SWEEP?

<b>Syntax</b>	SENSe:SWEEP?
<b>Description</b>	Queries the sweeping status of the laser
<b>Response</b>	ON if the laser is sweeping, OFF if the laser is stopped
<b>Example</b>	SENS:SWEEP?

#### 4.2.3.36 SENSE:LASer <C+L / O>

<b>Syntax</b>	SENSe:LASer <C+L   O>
<b>Description</b>	Selects the internal tuneable laser to operate with through the TLS OUT port. User can choose C+L or O band laser source. (Only available for BOSA O+C+L)
<b>Example</b>	SENS:LAS C+L

#### 4.2.3.37 SENSE:LASer ?

<b>Syntax</b>	SENSe:LASer ?
<b>Description</b>	Queries the actual internal laser that is being externally used for user purposes. (Only available for BOSA O+C+L)

<b>Comments</b>	This command is only available for the Tuneable Laser application.
<b>Response</b>	'Active internal tunable laser: C+L' or 'O'
<b>Example</b>	SENS:LAS ?

#### 4.2.3.38 SENSE:BAND <C+L / O / O+C+L>

<b>Syntax</b>	SENSe:BAND <C+L O O+C+L>
<b>Description</b>	Changes the measurement band to C+L, O or O+C+L band. BOSA application must be off to in order to change this option. Component Analyzer application does not support O+C+L option. (Only available for BOSA O+C+L)
<b>Example</b>	SENS:BAND O+C+L

#### 4.2.3.39 SENSE:BAND ?

<b>Syntax</b>	SENSe:BAND ?
<b>Description</b>	Queries the actual measured band. (Only available for BOSA O+C+L)
<b>Response</b>	'Active Measured Band: C+L' or 'O'
<b>Example</b>	SENS:BAND ?

### 4.2.4 INPut Subsystem Commands

#### 4.2.4.1 INPut:SPARameters

<b>Syntax</b>	INPut:SPARameters [IL RL IL&RL]
<b>Description</b>	Selects the measurement to display. The options are Insertion Loss (IL) and Return Loss (RL).
<b>Comments</b>	This command is only available for the Component Analyzer application
<b>Example</b>	INP:SPAR IL&RL

#### 4.2.4.2 INPut:SPARameters?

<b>Syntax</b>	INPut:SPARameters?
<b>Description</b>	Queries the type of measurement being displayed
<b>Response</b>	Possible responses are RL (return loss), IL (insertion loss) or IL&RL (insertion loss and return loss)

<b>Comments</b>	This command is only available for the Component Analyzer application
<b>Example</b>	INP:SPAR?

#### 4.2.4.3 INPut:POLarization

<b>Syntax</b>	<b>INPut:POLarization 1+2 1 2 1&amp;2</b>
<b>Description</b>	Selects the polarization measurement to be displayed.
<b>Argument</b>	It is possible to choose between the two orthogonal polarization states and display them separately or simultaneously.
<b>Example</b>	INP:POL 2

#### 4.2.4.4 INPut:POLarization

<b>Syntax</b>	<b>INPut:POLarization PDL MAX MIN SIMUL</b>
<b>Description</b>	Selects the polarization measurement to be displayed.
<b>Comments</b>	This command is only available for the Component Analyzer application with 230 option
<b>Example</b>	INP:POL PDL

#### 4.2.4.5 INPut:POLarization

<b>Syntax</b>	<b>INPut:POLarization INDEP 1 2 SIMUL</b>
<b>Description</b>	Selects the polarization measurement to be displayed.
<b>Comments</b>	This command is only available for the Component Analyzer application without 230 option
<b>Example</b>	INP:POL SIMUL

#### 4.2.4.6 INPut:POLarization?

<b>Syntax</b>	<b>INPut:POLarization?</b>
<b>Description</b>	Queries the current polarization configuration that is being displayed
<b>Response</b>	Depends on the current application and the hardware option installed.
<b>Example</b>	INP:POL?

#### 4.2.4.7 INPut:POLarization:MUELLermode

<b>Syntax</b>	<b>INPut:POLarization:MUELLermode &lt;ON OFF 1 0&gt;</b>
---------------	--

<b>Description</b>	Sets the advance polarization analysis on or off. Refer to the User Guide for further explaining.
<b>Comments</b>	This command is only available for the Component Analyzer application
<b>Example</b>	INP:POL:MUELL 0

#### 4.2.4.8 INPut:POLarization:MUELLermode?

<b>Syntax</b>	INPut:POLarization:MUELLermode?
<b>Description</b>	Queries the state of the polarization analysis. Refer to the User Guide for further explaining.
<b>Response</b>	It returns '1' or 'ON' when is set to the advance polarization mode and '0' or 'OFF' otherwise.
<b>Example</b>	INP:POL:MUELL?

#### 4.2.4.9 INPut:POWer?

<b>Syntax</b>	INPut:POWer?
<b>Description</b>	Queries total optical power being input to the BOSA. It can be also used in Option 210 for a readout of the Aux Input photodetector.
<b>Example</b>	INP:POW?

### 4.2.5 CALCulate Subsystem Commands

#### 4.2.5.1 CALCulate:MARKer:AOff

<b>Syntax</b>	CALCulate:MARKer:AOff
<b>Description</b>	Disable the marker.
<b>Example</b>	CALC:MARK:AOff
<b>Comments</b>	Referred to the active trace.

#### 4.2.5.2 CALCulate:MARKer:STaTe

<b>Syntax</b>	CALCulate:MARKer:STaTe ON OFF 1 0
<b>Description</b>	Enables or disables the marker.
<b>Argument</b>	'on' or '1' to turn the marker on, or 'off' or '0' to turn it off.
<b>Example</b>	CALC:MARK:STAT ON
<b>Comments</b>	Referred to the active trace.



#### 4.2.5.3 *CALCulate:MARKer:STATe?*

<b>Syntax</b>	<b>CALCulate:MARKer:STATe?</b>
<b>Description</b>	Queries the state of the marker.
<b>Response</b>	The current state of the marker: 'ON' or 'OFF'.
<b>Example</b>	CALC:MARK:STAT?
<b>Comments</b>	Referred to the active trace.

#### 4.2.5.4 *CALCulate:MARKer:MODE*

<b>Syntax</b>	<b>CALCulate:MARKer:MODE TRCK FIXX FIXXY</b>
<b>Description</b>	Sets the behaviour of the marker (and its offset if enabled).
<b>Argument</b>	The desired behaviour of the marker: 'TRCK', 'FIXX' or 'FIXXY'.
<b>Example</b>	CALC:MARK:MOD FIXX
<b>Comments</b>	Referred to the active trace.

#### 4.2.5.5 *CALCulate:MARKer:MODE?*

<b>Syntax</b>	<b>CALCulate:MARKer:MODE?</b>
<b>Description</b>	Queries the behaviour of the marker.
<b>Response</b>	The current behaviour of the marker: 'TRCK', 'FIXX' or 'FIXXY'.
<b>Example</b>	CALC:MARK:MOD?
<b>Comments</b>	Referred to the active trace.

#### 4.2.5.6 *CALCulate:MARKer:MAXimum*

<b>Syntax</b>	<b>CALCulate:MARKer:MAXimum</b>
<b>Description</b>	Moves the specified marker to the maximum power level.
<b>Example</b>	CALC:MARK:MAX
<b>Comments</b>	Referred to the active trace.

#### 4.2.5.7 *CALCulate:MARKer:MAXimum:NEXT*

<b>Syntax</b>	<b>CALCulate:MARKer:MAXimum:NEXT</b>
<b>Description</b>	Moves the marker to the next power maximum.
<b>Example</b>	CALC:MARK:MAX:NEXT
<b>Comments</b>	Referred to the active trace.

#### 4.2.5.8 *CALCulate:MARKer:MAXimum:RIGHT*

<b>Syntax</b>	<b>CALCulate:MARKer:MAXimum:RIGHT</b>
<b>Description</b>	Moves the marker to the next power maximum to the right.
<b>Example</b>	CALC:MARK:MAX:RIGHT
<b>Comments</b>	Referred to the active trace.

#### 4.2.5.9 *CALCulate:MARKer:MAXimum:LEFT*

<b>Syntax</b>	<b>CALCulate:MARKer:MAXimum:LEFT</b>
<b>Description</b>	Moves the marker to the next power maximum to the left.
<b>Example</b>	CALC:MARK:MAX:LEFT
<b>Comments</b>	Referred to the active trace.

#### 4.2.5.10 *CALCulate:MARKer:SCEN*

<b>Syntax</b>	<b>CALCulate:MARKer:SCEN</b>
<b>Description</b>	Moves the marker to the center of the screen.
<b>Example</b>	CALC:MARK:SCEN
<b>Comments</b>	Referred to the active trace.

#### 4.2.5.11 *CALCulate:MARKer:X*

<b>Syntax</b>	<b>CALCulate:MARKer:X &lt;numeric_value&gt; [NM PM GHZ THZ]</b>
<b>Description</b>	Moves the marker to the specified point.
<b>Argument</b>	'numeric_value' is the horizontal coordinate to put the marker on.
<b>Example</b>	CALC:MARK:X 1550 NM
<b>Comments</b>	Referred to the active trace.

#### 4.2.5.12 *CALCulate:MARKer:X?*

<b>Syntax</b>	<b>CALCulate:MARKer:X?</b>
<b>Description</b>	Queries the horizontal coordinate of the marker.
<b>Response</b>	The horizontal coordinate of the specified marker, expressed in units set by the READOUT function.
<b>Example</b>	CALC:MARK:X?
<b>Comments</b>	Referred to the active trace.

#### 4.2.5.13 CALCulate:MARKer:Y

<b>Syntax</b>	<b>CALCulate:MARKer:Y &lt;numeric_value&gt; [DBM MW]</b>
<b>Description</b>	Moves the marker to the specified point on the vertical axis.
<b>Argument</b>	'numeric_value' is the vertical coordinate to put the marker on.
<b>Example</b>	CALC:MARK:Y -20 DBM
<b>Comments</b>	Referred to the active trace.

#### 4.2.5.14 CALCulate:MARKer:Y?

<b>Syntax</b>	<b>CALCulate:MARKer:Y?</b>
<b>Description</b>	Queries the vertical coordinate of the marker.
<b>Response</b>	The vertical coordinate of the marker, expressed in the units set by the READOUT function.
<b>Example</b>	CALC:MARK:Y?
<b>Comments</b>	Referred to the active trace.

#### 4.2.5.15 CALCulate:MARKer:THREshold

<b>Syntax</b>	<b>CALCulate:MARKer:THREshold &lt;numeric_value&gt; [DB]</b>
<b>Description</b>	Sets the minimum power level for the search of peaks in peak search function
<b>Argument</b>	'numeric_value' is the depth of the peak search configuration.
<b>Example</b>	CALC:MARK:THRES-3 DB

#### 4.2.5.16 CALCulate:MARKer:THREshold?

<b>Syntax</b>	<b>CALCulate:MARKer:THREshold?</b>
<b>Description</b>	Queries the current peak excursion value for the peak search functions.
<b>Response</b>	The value for the peak threshold.
<b>Example</b>	CALC:MARK:THRES?

#### 4.2.5.17 CALCulate:MARKer:READout

<b>Syntax</b>	<b>CALCulate:MARKer:READout FREQuency WAVelenght</b>
<b>Description</b>	Sets the horizontal units of the marker
<b>Example</b>	CALC:MARK:READ FREQ

#### 4.2.5.18 CALCulate:MARKer:READout?

<b>Syntax</b>	<b>CALCulate:MARKer:READout?</b>
<b>Description</b>	Queries the horizontal units of the marker
<b>Example</b>	CALC:MARK:READ?

#### 4.2.5.19 CALCulate:MARKer:SRLevel

<b>Syntax</b>	<b>CALCulate:MARKer:SRLevel</b>
<b>Description</b>	Moves the reference level to the vertical value of the marker.
<b>Example</b>	CALC:MARK:SRL
<b>Comments</b>	Referred to the active trace. Activate marker on desired trace first, otherwise the command will return error.

#### 4.2.5.20 CALCulate:MARKer:POLarization?

<b>Syntax</b>	<b>CALCulate:MARKer:FUNCTION:DELTA:POLarization?</b>
<b>Description</b>	In polarimetry extension, queries the S parameters (S1, S2, S3) of the marker in ASCII format.
<b>Example</b>	CALC:MARK:POL?
<b>Comments</b>	Values separated by commas.

#### 4.2.5.21 CALCulate:MARKer:FUNCTION:DELTA

<b>Syntax</b>	<b>CALCulate:MARKer:FUNCTION:DELTA [:STATE] ON OFF 1 0</b>
<b>Description</b>	Turns the delta marker function on or off.
<b>Example</b>	CALC:MARK:FUNC:DELTA ON
<b>Comments</b>	Referred to the active trace.

#### 4.2.5.22 CALCulate:MARKer:FUNCTION:DELTA?

<b>Syntax</b>	<b>CALCulate:MARKer:FUNCTION:DELTA [:STATE]?</b>
<b>Description</b>	Queries the state of the delta marker function.
<b>Example</b>	CALC:MARK:FUNC:DELT?
<b>Comments</b>	Referred to the active trace.

#### 4.2.5.23 CALCulate:MARKer:FUNCTION:DELTA:RESet

<b>Syntax</b>	<b>CALCulate:MARKer:FUNCTION:DELTA:RESet</b>
<b>Description</b>	Sets the marker offset to the current position of the marker.

<b>Example</b>	CALC:MARK:MARK:FUNC:DELT:RES
<b>Comments</b>	Referred to the active trace.

#### 4.2.5.24 CALCulate:MARKer:FUNCTION:DELTa:X:OFFSet?

<b>Syntax</b>	CALCulate:MARKer:FUNCTION:DELTa:X:OFFSet?
<b>Description</b>	Queries the horizontal offset between the marker and its reference.
<b>Example</b>	CALC:MARK:FUNC:DELT:X:OFFS?
<b>Comments</b>	Referred to the active trace.

#### 4.2.5.25 CALCulate:MARKer:FUNCTION:DELTa:X:REFerence?

<b>Syntax</b>	CALCulate:MARKer:FUNCTION:DELTa:X:REFerence?
<b>Description</b>	Queries the horizontal value of the reference marker.
<b>Example</b>	CALC:MARK:FUNC:DELT:X:REF?
<b>Comments</b>	Referred to the active trace.

#### 4.2.5.26 CALCulate:MARKer:FUNCTION:DELTa:Y:OFFSet?

<b>Syntax</b>	CALCulate:MARKer:FUNCTION:DELTa:Y:OFFSet?
<b>Description</b>	Queries the vertical offset between the marker and its reference.
<b>Example</b>	CALC:MARK:FUNC:DELT:Y:OFFS?
<b>Comments</b>	Referred to the active trace.

#### 4.2.5.27 CALCulate:MARKer:FUNCTION:DELTa:Y:REFerence?

<b>Syntax</b>	CALCulate:MARKer:FUNCTION:DELTa:Y:REFerence?
<b>Description</b>	Queries the vertical value of the reference marker.
<b>Example</b>	CALC:MARK:FUNC:DELT:Y:REF?
<b>Comments</b>	Referred to the active trace.

#### 4.2.5.28 CALCulate:MARKer:FUNCTION:DELTa:POL?

<b>Syntax</b>	CALCulate:MARKer:FUNCTION:DELTa:POLarization?
<b>Description</b>	In polarimetry extension, queries the S parameters (S1, S2, S3) of the offset in ASCII format.
<b>Example</b>	CALC:MARK:FUNC:DELT:POL?
<b>Comments</b>	Values separated by commas.

#### 4.2.5.29 CALCulate:MARKer:FUNCTION:DELTA:ANG?

<b>Syntax</b>	<b>CALCulate:MARKer:FUNCTION:DELTA:ANGLE?</b>
<b>Description</b>	In polarimetry extension, queries the angle between marker and offset in ASCII format.
<b>Example</b>	CALC:MARK:FUNC:DELT:ANG?
<b>Comments</b>	

#### 4.2.5.30 CALCulate:MAXimum[:STATe]

<b>Syntax</b>	<b>CALCulate:MAXimum[:STATe] 0 1 ON OFF</b>
<b>Description</b>	Enables or disables the Max Hold function.
<b>Argument</b>	'1' or 'ON' to enable the Max Hold function; '0' or 'OFF' to disable it.
<b>Example</b>	CALC:MAX ON

#### 4.2.5.31 CALCulate:MAXimum[:STATe]?

<b>Syntax</b>	<b>CALCulate:MAXimum[:STATe]?</b>
<b>Description</b>	Queries the Max Hold function state.
<b>Response</b>	'ON' if it is enabled, 'OFF' if not.
<b>Example</b>	CALC:MAX?

#### 4.2.5.32 CALCulate:MINimum[:STATe]

<b>Syntax</b>	<b>CALCulate:MINimum[:STATe] 0 1 ON OFF</b>
<b>Description</b>	Enables or disables the Min Hold function.
<b>Argument</b>	'1' or 'ON' to enable the Min Hold function; '0' or 'OFF' to disable it.
<b>Example</b>	CALC:MIN ON

#### 4.2.5.33 CALCulate:MINimum[:STATe]?

<b>Syntax</b>	<b>CALCulate:MINimum[:STATe]?</b>
<b>Description</b>	Queries the Min Hold function state.
<b>Response</b>	'ON' if it is enabled, 'OFF' if not.
<b>Example</b>	CALC:MIN?

#### 4.2.5.34 CALCulate:TPOWer ON/OFF

Syntax	CALCulate:TPOWer 0 1 ON OFF
Description	Enables or disables the power integral function.
Argument	'1' or 'ON' to enable the Power Integral function; '0' or 'OFF' to disable it.
Example	CALC:TPOW ON

#### 4.2.5.35 CALCulate:TPOWer:[DATA]?

Syntax	CALCulate:TPOWer:[DATA]?
Description	Returns the power integral in mW or dBm.
Example	CALC:TPOW?

#### 4.2.5.36 CALCulate:TPOWer:IRANge:UPPer

Syntax	CALCulate:TPOWer:IRANge:UPPer<numeric_value> [nm pm THz GHZ]
Description	Sets the upper bound for power integration.
Argument	'numeric_value' is the value for the upper bound. Units are optional.
Example	CALC:TPOW:IRAN:UPP 1555 NM

#### 4.2.5.37 CALCulate:TPOWer:IRANge:UPPer?

Syntax	CALCulate:TPOWer:IRANge:UPPer?
Description	Returns the upper bound for power integration in nm or GHz.
Example	CALC:TPOW:IRAN:UPP?

#### 4.2.5.38 CALCulate:TPOWer:IRANge:LOWer

Syntax	CALCulate:TPOWer:IRANge:LOWer [nm pm THz GHZ] <numeric_value>
Description	Sets the lower bound for power integration.
Argument	'numeric_value' is the value for the lower bound. Units are optional.
Example	CALC:TPOW:IRAN:LOW 1550 NM

Syntax	CALCulate:TPOWer:IRANge:LOWer?
--------	--------------------------------

<b>Description</b>	Returns the lower bound for power integration in nm or GHz.
<b>Example</b>	CALC:TPOW:IRAN:LOW?

#### 4.2.5.39 CALCulate:TPOWer:IRANge:LOWer?

#### 4.2.5.40 CALCulate:AUXINput:POWER?

<b>Syntax</b>	CALCulate:TPOWer:IRANge:LOWer?
<b>Description</b>	Returns the power meter measurement from the Aux Input port in dBm. Only available in the tunable laser application.
<b>Example</b>	CALC:AUXIN:POW?

## 4.2.6 TRACe Subsystem Commands

#### 4.2.6.1 TRACe[:DATA]:COUNT?

<b>Syntax</b>	TRACe[:DATA]:COUNT?
<b>Description</b>	Queries the number of points displayed on the screen.
<b>Response</b>	The number of points showed in the BOSA screen. Each point represent a 64 bit numeric value.
<b>Example</b>	TRAC:COUNT?

#### 4.2.6.2 TRACE[:DATA]?

<b>Syntax</b>	TRACe[:DATA]?
<b>Description</b>	Queries the whole spectrum measured by the BOSA.
<b>Response</b>	The entire trace displayed in the BOSA screen. Use FORMAT[:DATA] to specify the format of the returned data. When format is set to REAL points are 64 bit in length. Use TRACE:COUNT? to know the number of points returned. When format is set to ASCII the values are separated by commas, the whole returned string finishes with "\r\n".



<b>Comments</b>	This command only works when the BOSA is controlled remotely.
<b>Example</b>	TRAC?

#### 4.2.6.3 TRACe[:DATA]:MAXimum:X?

<b>Syntax</b>	TRACe[:DATA]:MAXimum:X?
<b>Description</b>	Queries the X coordinate of the absolute maximum.
<b>Response</b>	The X coordinate of the absolute maximum.
<b>Example</b>	TRAC:MAX:X?

#### 4.2.6.4 TRACe[:DATA]:MAXimum:Y?

<b>Syntax</b>	TRACe[:DATA]:MAXimum:Y?
<b>Description</b>	Queries the Y coordinate of the absolute maximum.
<b>Response</b>	The Y coordinate of the absolute maximum.
<b>Example</b>	TRAC:MAX:Y?

### 4.2.7 FORMat Subsystem Commands

#### 4.2.7.1 FORMat[:DATA]

<b>Syntax</b>	FORMat[:DATA] <ASCII REAL>,<length>
<b>Description</b>	Sets the output data format for TRACe[:DATA][:Y]?. Available formats are either to ASCII with a given length or REAL (64 bits).
<b>Argument</b>	'ASCII' or 'REAL' for the format and the 'length' of the string in case of the ASCII option.
<b>Example</b>	FORM ASCII,5

#### 4.2.7.2 FORMat[:DATA]?

<b>Syntax</b>	FORMat[:DATA]?
<b>Description</b>	Returns the current output data format for TRACe[:DATA][:Y]?
<b>Response</b>	Example: 'ASCII,6' or 'REAL,6'
<b>Example</b>	FORM?

## 4.2.8 MMEMory Subsystem Commands

### 4.2.8.1 MMEMory:STORe:TRACe

<b>Syntax</b>	<b>MMEMory:STORe:TRACe &lt;name&gt;.&lt;type&gt;,[ink_saving_mode]</b>
<b>Description</b>	Stores the specified trace in a file in the specified format.
<b>Arguments:</b>	<p>'name' is the path of the new file (without the extension) and 'type' can be: bdf, txt, csv, jpg, bmp, gif or tif. Finally, 'ink_saving_mode' is an optional Boolean (1, 0, ON or OFF) parameter only for image formats. If it is true the image will be saved with ink saving mode colors. If it is omitted, the image will display the original BOSA colors set up.</p> <p>NOTE: the root directory is 'User Files'. The file name is relative to the root directory.</p>
<b>Example</b>	MMEM:STOR:TRAC myImages\image1.jpg,1

### 4.2.8.2 MMEMory:DELeTe:TRACe

<b>Syntax</b>	<b>MMEMory:DELeTe:TRACe &lt;name&gt;</b>
<b>Description</b>	Deletes the specified file
<b>Arguments:</b>	'name' is the file name to be deleted. The root directory is 'User Files'. The file name is relative to the root directory.
<b>Example</b>	MMEM:DEL:TRAC gpr.csv

### 4.2.8.3 MMEMory:LOAD:TRACe

<b>Syntax</b>	<b>MMEMory:LOAD:TRACe &lt;M1 M2 M3 M4&gt;,&lt;name&gt;</b>
<b>Description</b>	Loads the specified file
<b>Arguments:</b>	'M1' to 'M4' sets the trace in which the BOSA software will load the file, 'name' is the file name to be loaded. The root directory is 'User Files'. The file name is relative to the root directory.
<b>Example</b>	MMEM:LOAD:TRAC M2,gpr.csv





# 5 Programming examples

## 5.1 Ethernet Interface

The following example of C code shows how to communicate with the BOSA using the Ethernet Interface. A query is sent to the BOSA and its response is read.

```
// SocketsCExample.cpp : Uses the ethernet communication
// open a socket, send *IDN? query and print the response
#include "stdafx.h"

// include socket library
#include <winsock2.h>
#pragma comment(lib, "ws2_32.lib")

#define LEN_REP 255 //length of the response buffer

int main(int argc, char* argv[])
{
    WSADATA WSAData;
    SOCKET sock;
    SOCKADDR_IN sin;
    char response[LEN_REP];
    int byte_read;
    int i;

    //initialize WINSOCK
    WSStartup(MAKEWORD(2,0), &WSAData);

    //Take the address of the instrument
    sin.sin_addr.s_addr = inet_addr("192.168.127.79");
    sin.sin_family = AF_INET;
    sin.sin_port = htons(10000);

    if((sock=socket(sin.sin_family,SOCK_STREAM,0)) < 0)
    {
        perror("Creation socket failed");
        return(1);
    }

    if(connect(sock, (SOCKADDR *)&sin, sizeof(sin))<0)
    {
        perror("Connect failed");
        return(1);
    }
}
```

```
//send the command
send(sock, "*IDN?", strlen("*IDN?"), 0);

//-1: leave space for null terminator

if((byte_read=recv(sock, response, LEN_REP-1, 0))<=0)
{
    perror("Receive error");
    return(1);
}
response[byte_read] = 0x00; // terminate the string

printf("%s",response);
closesocket(sock);
WSACleanup();

return 0;
```

## 5.2 Macro Editor Tool

### 5.2.1 Performing auto measurements of several peaks using markers

The following example shows how to use markers to perform automatic searching and measurements of several peaks in the spectrum.

```
////////////////////////////////////
//          MULTIPEAK MEASUREMENT          //
////////////////////////////////////

//      DESCRIPTION

//      This macro will search for peaks present in the configured span,
//      save their wavelength and power in an array and finally store the
//      array values in a file.

//      DECLARATION AND INITIALIZATION OF VARIABLES

CENTERWL = 1550                                // The center wavelength in nm
SPAN = 5                                         // The desired span, in nm
REFLEVEL = -50                                  // The Reference level
DATE = FORMATDATE(NOW()),"MMMM_dd hh_mm")      // The actual date
FILE = "PeaksMeasurement\results " & DATE & ".txt" // Results will be saved in this file
PEAKPOSITION = ARRAY(20)                       // It will contain the peak positions
PEAKPOWER = ARRAY(20)                          // It will contain the peak powers

COUNT=0
```

```
//      INSTRUCTIONS

//      Configuration of the BOSA: Center and span
SCPI("SENS:WAV:CENT" & " " & CENTERWL & " " & "NM")
SCPI("SENS:WAV:SPAN" & " " & SPAN & " " & "NM")

//      Configuration of the BOSA: Peak Search parameters
//      Markers will only find peaks with a peak excursion greater than
//      'pexcursion' at a distance of 'pwidth' from the wavelength of the
//      maximum of the peak
SCPI("CALC:MARK1:PEXCURSION 10 DB")
SCPI("CALC:MARK1:PWIDITH 0.0002 NM")

//      Configuration of the BOSA: Reference Level.
//      The reference level is used as threshold for peak search.
//      Peaks with amplitude lower than refllevel will not be found
SCPI("DISP:TRAC:Y:RLEV" & " " & REFLEVEL & " " & "DBM")

//      Activate the marker 1
SCPI("CALC:MARK1:STAT ON")
//      Move the marker 1 to the starting wavelength
SCPI("CALC:MARK1:X STAR")
//      Move the marker 1 to the first peak
SCPI("CALC:MARK1:MAX:NEXT")
WAIT(200)

//      Ask the BOSA the peak parameters
POSITION = TODOUBLE(SCPI("CALC:MARK1:X?"))
POWER = TODOUBLE(SCPI("CALC:MARK1:Y?"))

//      MEASUREMENT LOOP
REPEAT
    PEAKPOSITION[COUNT] = POSITION
    PEAKPOWER[COUNT] = POWER
    COUNT = COUNT+1
    SCPI("CALC:MARK1:MAX:NEXT")
    WAIT(200)
    POSITION = TODOUBLE(SCPI("CALC:MARK1:X?"))
    POWER = TODOUBLE(SCPI("CALC:MARK1:Y?"))
    WHILE(POSITION<PEAKPOSITION[COUNT-1]&&COUNT<LENGTH(PEAKPOSITION))

    TRIM(PEAKPOSITION)
    TRIM (PEAKPOWER)

//      SAVING DATA TO A FILE
NUMFORMAT = "0.0000"
FWRITE(FILE,"SAMPLE MACRO FOR MULTIPLE PEAKS AUTOMATED MEASUREMENT")
FWRITE(FILE,"")
FWRITE(FILE,"PEAK#  WL(NM)  POW(DBM)")
j=0
REPEAT
    FWRITE(FILE,TOSTRING(J+1))
    FWRITE(FILE," " & FORMATNUM(PEAKPOSITION[J] , NUMFORMAT ))
    FWRITE(FILE," " & FORMATNUM(PEAKPOWER[J] , NUMFORMAT ))
    J=J+1
WHILE(J<LENGTH("PEAKPOWER"))
```

## 5.2.2 Wavelength shift measurement on a CW laser

The following example shows how to perform an automatic measurement of the wavelength shift on a CW laser over time.

```

////////////////////////////////////
//  WAVELENGTH SHIFT MEASUREMENT ON A CW LASER  //
////////////////////////////////////

//  DESCRIPTION

//  The following program measures the wavelength emission of a CW laser
//  over time. A wavelength measurement is taken and stored in a data
//  array at the specified time interval, finally the results are plot in a chart.

//  DECLARATION AND INITIALIZATION OF VARIABLES

NUMMEAS=100           // Number of measurements to be taken
TIMEINTERVAL=10       // Time interval between measurements in seconds
INDEX=0

DATA=ARRAY(NUMMEAS)

//  INSTRUCTIONS

//  An auto-measurement is performed to locate the output laser wavelength
SCPI("DISP:TRAC:X:AUTO")
WAIT(20000)

//  MEASUREMENT LOOP

REPEAT
//  It queries the wavelength of the power maximum
WMAX=SCPI("TRACE:DATA:MAXIMUM:X?")
//  The wavelength is store in the data array
DATA[INDEX]=TODouble(WMAX)
//  The measurement center wavelength is set to the wavelength of the
//  measured maximum. It avoids the laser to shift out the measurement span.
SCPI("SENS:WAV:CENTER " & WMAX & " NM")
//  The time interval between measurements is used in milliseconds
WAIT(TIMEINTERVAL*1000)
INDEX=INDEX+1
WHILE (INDEX<NUMMEAS)

//  The DATA array is plot in a chart
PLOT(DATA)

SORT(DATA)
MAXIMUM=DATA[NUMMEAS-1]
MINIMUM=DATA[0]
DELTA=MAXIMUM-MINIMUM

SCPI("INST:STAT OFF")

//  It is possible to reach the numeric data by following the path File -> Save as Text on the
//  chart. A txt file is
//  saved with all the numeric values plotted on the chart.

```



### 5.2.3 Controlling an external device via GPIB

The following example shows how to manage an external device through GPIB using its own SCPI commands. Note that these commands depend on each device and they could be different for each model.

```

////////////////////////////////////
//          MANAGING AN EXTERNAL TLS VIA GPIB          //
////////////////////////////////////

//      DESCRIPTION

//      The following program measures the wavelength and the power of
//      an external TLS which is set to perform stepped sweeps.

//      DECLARATION AND INITIALIZATION OF VARIABLES

WSTART= 1540
WSTOP = 1560
INTERVAL = 2
POWER = 0
WLENGTH = 0
COUNT = 0
LAMBDA = 0
FILE = "externalTLS.txt"           // Results will be saved in this file
IDTLS = 10                        // The GPIB ID number of the TLS
SLOTTLS = "0"                    // The slot number of the TLS???

WAVELENGTHTHEO = ARRAY(11)        // 11 = size = 1 + (wstop-wstart)/interval
WAVELENGTHREAL = ARRAY(11)
POWERREAL = ARRAY(11)

//      INSTRUCTIONS

//      Configuration of the BOSA: Start and stop
SCPI("SENS:WAV:START" & " " & (WSTART-INTERVAL)& " " & "NM")
SCPI("SENS:WAV:STOP" & " " & (WSTOP+INTERVAL) & " " & "NM")

//      The BOSA is configured to measure at the highest resolution
SCPI("SENS:WAV:RES HIGH")

LAMBDA = WSTART

//      Configuration of the external TLS: Lambda and power
SENDGPIB(IDTLS, "SOUR"&SLOTTLS&":WAV "& FORMATNUM(LAMBDA,"0.###")&" NM")
SENDGPIB(IDTLS, "SOUR"&SLOTTLS&":POW "& FORMATNUM(POWER,"0.0")&" DBM")
//      Turn on the external TLS
SENDGPIB(IDTLS, "SOUR"&SLOTTLS&":POW:STAT 1")

WAIT(7500)

//      It queries the wavelength of the maximum
WLENGTH = TODOUBLE(SCPI("TRACE:DATA:MAXIMUM:X?"))
//      It queries the power of the maximum
POWER = TODOUBLE(SCPI("TRACE:DATA:MAXIMUM:Y?"))

//      MEASUREMENT LOOP

REPEAT
    WAVELENGTHTHEO[COUNT] = LAMBDA
    WAVELENGTHREAL[COUNT] = WLENGTH
    POWERREAL[COUNT] = POWER
    COUNT = COUNT+1

    LAMBDA = LAMBDA +INTERVAL

```

```

//      Configuring the external TLS: Lambda
SENDGPB(IDTLS, "SOUR"&SLOTTLS&":WAV "& FORMATNUM(LAMBDA,"0.###")&" NM")

WAIT(5000)

//      It queries the wavelength of the maximum
WLENGTH = TODOUBLE(SCPI("TRACE:DATA:MAXIMUM:X?"))
//      It queries the the power of the maximum
POWER = TODOUBLE(SCPI("TRACE:DATA:MAXIMUM:Y?"))
WHILE (LAMBDA<=WSTOP)

//      Turn off the external TLS
SENDGPB(IDTLS, "SOUR"&SLOTTLS&":POW:STAT 0")

//      SAVING DATA TO A FILE

NUMFORMAT = "0.0000"
FWRITE(FILE,"SAMPLE MACRO FOR MANAGING EXTERNAL TLS")
FWRITE(FILE,"")
FWRITE(FILE,"WL THEO (NM)   WL REAL (NM)   POWER(DBM)")
J=0
REPEAT
    FWRITE(FILE,"   "& FORMATNUM(WAVELENGTHTHEO[J] , NUMFORMAT ))
    FWRITE(FILE,"   "& FORMATNUM(WAVELENGTHREAL[J] , NUMFORMAT ))
    FWRTTELN(FILE," "& FORMATNUM(POWERREAL[J] , NUMFORMAT ))
    J=J+1
WHILE(J<LENGTH(POWERREAL ))

```