

Database connection type: PostgreSQL
Host name: sysbio-db.kybi
Database: experiments_v2
Schema: public
User: markov

Open new experiment

Change database


database connection

date	time	type	info	hardware
2018.03.13	12:22	VO2 ADP t.	channel-3	strathkelvin
2018.03.13	12:22	VO2 ADP t.	channel-6	strathkelvin
2018.03.13	13:53	VO2 ADP t.	channel-3	strathkelvin
2018.03.16	12:23	VO2 ADP t.	channel-6	strathkelvin
2018.03.16	14:04	VO2 ADP t.	channel-2	strathkelvin
2018.03.22	12:54	VO2 ADP t.	channel-2	strathkelvin
2018.03.22	12:54	VO2 ADP t.	channel-3	strathkelvin
2018.03.22	12:54	VO2 ADP t.	channel-6	strathkelvin
2018.03.23	12:38	VO2 ADP t.	channel-2	strathkelvin
2018.03.23	12:38	VO2 ADP t.	channel-3	strathkelvin
2018.03.23	12:38	VO2 ADP t.	channel-6	strathkelvin
2018.04.12	12:39	VO2 ADP t.	channel-3	strathkelvin
2018.04.12	12:39	VO2 ADP t.	channel-6	strathkelvin
2018.04.16	12:04	VO2 ADP t.	channel-3	strathkelvin
2018.04.16	12:04	VO2 ADP t.	channel-6	strathkelvin
2018.04.17	12:14	VO2 ADP t.	channel-3	strathkelvin
2018.04.17	12:14	VO2 ADP t.	channel-6	strathkelvin
2018.04.18	12:55	VO2 ADP t.	channel-3	strathkelvin
2018.04.18	12:55	VO2 ADP t.	channel-6	strathkelvin
2018.04.20	12:39	VO2 ADP t.	channel-3	strathkelvin
2018.04.20	12:39	VO2 ADP t.	channel-6	strathkelvin
2018.05.02	12:43	VO2 ADP t.	channel-3	strathkelvin
2018.05.02	12:43	VO2 ADP t.	channel-6	strathkelvin
2018.05.03	12:30	VO2 ADP t.	channel-3	strathkelvin
2018.05.03	12:30	VO2 ADP t.	channel-6	strathkelvin
2018.05.04	11:39	VO2 ADP t.	channel-3	strathkelvin
2018.05.04	11:39	VO2 ADP t.	channel-6	strathkelvin
2018.05.08	11:52	VO2 ADP t.	channel-3	strathkelvin
2018.05.08	11:52	VO2 ADP t.	channel-6	strathkelvin
2018.05.09	11:59	VO2 ADP t.	channel-2	strathkelvin
2018.05.29	11:59	VO2 ADP t.	channel-3	strathkelvin
2018.05.29	11:59	VO2 ADP t.	channel-6	strathkelvin
2018.05.30	11:51	VO2 ADP t.	channel-3	strathkelvin
2018.05.30	11:51	VO2 ADP t.	channel-6	strathkelvin
2018.06.01	12:22	VO2 ADP t.	channel-3	strathkelvin
2018.06.01	12:22	VO2 ADP t.	channel-6	strathkelvin
2018.06.04	12:21	VO2 ADP t.	channel-3	strathkelvin
2018.06.04	12:21	VO2 ADP t.	channel-6	strathkelvin
2018.06.05	12:06	VO2 ADP t.	channel-3	strathkelvin
2018.06.05	12:06	VO2 ADP t.	channel-6	strathkelvin
2018.06.05	14:11	VO2 ADP t.	channel-3	strathkelvin
2018.06.05	14:11	VO2 ADP t.	channel-6	strathkelvin
2018.06.18	14:52	VO2 ADP t.	channel-2	strathkelvin
2018.06.18	14:52	VO2 ADP t.	channel-3	strathkelvin
2018.06.19	13:52	VO2 ADP t.	channel-1	strathkelvin
2018.06.19	13:52	VO2 ADP t.	channel-4	strathkelvin
2018.06.20	14:16	VO2 ADP t.	channel-1	strathkelvin
2018.06.26	12:13	VO2 ADP t.	channel-4	strathkelvin
2018.06.28	13:58	VO2 ADP t.	channel-1	strathkelvin
2018.06.28	13:58	VO2 ADP t.	channel-4	strathkelvin

read/read-write mode switch

Read only mode. Note: When read-only, no changes will be saved


importing data



ROI list:


- 1
- 2
- 3
- 4
- 5
- 6
- 7
- 8
- 9
- 10

ID



Event name: Cyt c

parameter	value	unit
Respiration rate	6.3314	umol/l/min
p	< 10^-12	
r-squared	0.99978	
r	-0.99989	
std err	0.005686	



parameter	value	unit
ACR: (Vmax-V0)/V0	8.4635	
Flux control ratio: V0/Vmax	0.10567	
V0 (after subtracting Vbg)	0.62833	umol/l/min
Vbg (background rate)	0.38517	umol/l/min
V maximal recorded (after subtracting Vbg)	5.9462	umol/l/min
Km	0.36405	mM or other unit
Vmax	5.8738	umol/l/min